

KLS martin
GROUP



marCut[®]

The bipolar KLS Martin Scissors
for use in Open Surgery




marCut®

The Bipolar KLS Martin Scissors for Use in Open Surgery


The bipolar technique

When using the bipolar electro-surgical technique, the flow of current is limited to a defined tissue area. Only the tissue between the two blades (tips) of the scissors is part of the electric circuit.


This requires significantly less power than mono-polar applications. In other words, the bipolar technique means more patient safety!




Precise cutting and simultaneous hemostasis – in a single surgical step.



Coagulation is possible in any blade position. Thus, even extended superficial hemorrhages can be easily stopped.



The bipolar scissors can also be used for point coagulation (instead of forceps) to stop locally limited hemorrhages.




With the scissors opened, vessel coagulation is possible before performing the cut.

Safe and universally applicable Accessories for Electrosurgery

marCut®

The bipolar KLS Martin scissors for use in open surgery



Scissor blade with ceramic coating

To effectively protect the ceramic-coated blade of the marCut® scissors during instrument processing, it is recommended to use the protector (80-570-00-04).

Design of the bipolar KLS Martin scissors

The two blades of the scissors are electrically insulated from each other by a ceramic element. The specially hardened surface of the other blade guarantees a consistently sharp cutting edge.



Good reasons at a glance


- Tissue-friendly and safe treatment as an alternative to monopolar surgery
- Durable ceramic coating; the specially hardened surface of the other blade ensures a consistently sharp edge
- Rugged metal joint (pivot)
- Can be resharpened repeatedly
- Can be used across a wide range of indications
- Time savings due to safe coagulating and precise cutting in a single step
- Both the scissors and the cables can be steam-sterilized (up to 134°C / 2 bar)
- Existing bipolar connecting cables can be used as well
- The marCut® scissors can be used with any KLS Martin HF generator capable of producing bipolar currents; Compatible non-KLS-Martin generators on request
- Connection of the bipolar cable to only one shank (or handle) facilitates the surgeon's work

Bipolar Scissors Current for electrosurgically assisted mechanical Cutting

With its maxium®, the latest generation of HF units, KLS Martin offers you a special scissors current: Forfex. This current has been specially designed for use with our bipolar marCut® scissors. It supports the mechanical cut of the scissors with bipolar high-frequency current that enables simultaneous hemostasis.



Application: colon resection



Application: thyroidectomy



The maxium® generator's Forfex current has been specially adapted to the marCut® bipolar scissors.

Indication examples

General surgery

- Laparotomy
- Intestinal surgery
- Esophageal interventions
- Colorectal interventions
- Pancreatic interventions
- Hernias
- Adhesions

Gynecology

- Abdominal hysterectomy
- Vaginal hysterectomy
- Mamma dissection
- Axillary lymph node dissection
- Adhesions

Thoracic surgery

- Thoracotomy
- Pneumonectomy
- Adhesions

Plastic surgery

- Flap plasty
- Mamma reconstruction

Cardiac surgery

- Excisional biopsies
- Mammary dissection
- Saphenous vein preparation

ENT


- Tonsillectomies
- Neck dissections
- Thyroidectomies
- Oral cavity resection

Types of Scissors and Connecting Cables

Different sizes give you choice – just the right instrument for the task at hand. We recommend using the protector (80-570-00-04) for processing the scissors.

Users can choose between the standard version and two versions of the slim model – the marCut® SL (Slim Line) and the new marCut® SSL (Super Slim Line).

Standard version





max. 600 V_p

80-287-33-04 80-287-53-04
3 m/10 ft. 5 m/16 ft.
KLS Martin and Berchtold units
→ bipolar instruments



max. 600 V_p

80-286-30-04 80-286-50-04
3 m/10 ft. 5 m/16 ft.
maxium® "e"-Version/ICC and
ACC/VIO Erbe units
→ bipolar instruments




max. 600 V_p

80-287-89-04 80-287-90-04
3 m/10 ft. 5 m/16 ft.
maxium® "i"-Version/ME MB 3 "i"-Version,
ME MB 2 "i"-Version/Valleylab units
→ bipolar instruments



max. 600 V_p

80-287-72-04
5 m/16 ft.
Codman and Valleylab units
→ bipolar instruments



80-570-00-04

Protector for instrument processing.
We recommend using the protector for
processing the scissors.

Slim Line Versión estilizada



80-571-18-04
marCut® SL
18 cm/7"



80-571-21-04
marCut® SL
21 cm/8 1/4"



80-571-23-04
marCut® SL
23 cm/9"



80-571-28-04
marCut® SL
28 cm/11"

Types of Scissors and Connecting Cables

Different sizes give you choice – just the right instrument for the task at hand. We recommend using the protector (80-570-00-04) for processing the scissors.

Users can choose between the standard version and two versions of the slim model – the marCut® SL (Slim Line) and the new marCut® SSL (Super Slim Line).

Super Slim Line Super slim model





1/3 max. 600 V_p

80-287-33-04 80-287-53-04
 3 m/10 ft. 5 m/16 ft.
 KLS Martin and Berchtold units
 → bipolar instruments



1/3 max. 600 V_p

80-286-30-04 80-286-50-04
 3 m/10 ft. 5 m/16 ft.
 maxium® "e"-Version/ICC and
 ACC/VIO Erbe units
 → bipolar instruments




1/3 max. 600 V_p

80-287-89-04 80-287-90-04
 3 m/10 ft. 5 m/16 ft.
 maxium® "i"-Version/ME MB 3 "i"-Version,
 ME MB 2 "i"-Version/Valleylab units
 → bipolar instruments



1/3 max. 600 V_p

80-287-72-04
 5 m/16 ft.
 Codman and Valleylab units
 → bipolar instruments



80-570-00-04
 Protector for instrument processing.
 We recommend using the protector for
 processing the scissors.

Should you have any queries ...
 ... we'll be glad to be at your Service!



Our HF Accessories Catalog
 contains the complete range
 of KLS Martin HF instruments

You can reach us in person either via
 e-mail or through our customer hotline.
 E-mail: info@klsmartin.com
 Hotline: ++49-7461-706-0

KLS Martin Group

KLS Martin Australia Pty Ltd.

Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Medical (Shanghai) International Trading Co. Ltd.

Shanghai · China
Tel. +86 21 5820 6251
china@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
italia@klsmartin.com

Nippon Martin K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
nippon@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 505 7838
malaysia@klsmartin.com

KLS Martin de México S.A. de C.V.

Mexico City · Mexico
mexico@klsmartin.com

Martin Nederland/Marned B.V.

Huizen · The Netherlands
Tel. +31 35 523 45 38
nederland@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Moscow · Russia
Tel. +7 499 792-76-19
russia@klsmartin.com

KLS Martin Taiwan Ltd.

Taipei 106 · Taiwan
Tel. +886 2 2325 3169
taiwan@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.

London · United Kingdom
Tel. +44 1189 000 570
uk@klsmartin.com

KLS Martin LP

Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

Do you know how to get all important information
about the KLS Martin Energy Devices?

Please download KLS Martin App Energy Devices!
The App is available for Android and iOS.



[https://itunes.apple.com/de/app/
klsmartin-energy-devices/id1198171415?l=de&ls=1&mt=8](https://itunes.apple.com/de/app/klsmartin-energy-devices/id1198171415?l=de&ls=1&mt=8)



[https://play.google.com/store/apps/
details?id=com.klsmartin.energydevices](https://play.google.com/store/apps/details?id=com.klsmartin.energydevices)

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany
P.O. Box 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com

