

Care and Cleaning

Information for the BK Medical Product Range



LEGAL MANUFACTURER

BK Medical ApS

Mileparken 34

2730 Herlev

Denmark

Tel.:+45 4452 8100 / Fax:+45 4452 8199

www.bkmedical.com

Email: info@bkmedical.com

The serial number of a BK Medical product contains information about the year of manufacture.

New disinfection and sterilization methods are constantly being developed, and we work to make our products compatible with as many methods as possible. You can find the latest information about caring for our products, including disinfection and sterilization compatibility, on our website. New information may have been added since you received this book.

BK Medical Customer Satisfaction

Input from our customers helps us improve our products and services. Your opinions are important to us. You are always welcome to contact us via your BK Medical representative or by contacting us directly.

Should any serious incident occur in relation to this device, you should report this to the manufacturer and your local competent authority.










Contents








| | |
|--|-----------|
| Chapter 1: Warnings and Cautions | 5 |
| Warnings | 6 |
| Cautions | 7 |
| | 11 |
| Chapter 2: Markings on the Systems, Transducers and Accessories | 11 |
| Safety Symbols and Information on the Equipment | 12 |
| Chapter 3: Introduction & Safety | 17 |
| Introduction | 18 |
| Warnings, Cautions, Notes | 18 |
| General Safety | 18 |
| bkSpecto | 19 |
| The bk3000/bk5000/bk3500 System and UA2361 Remote Control | 19 |
| The Flex Focus 1202 System and UA1237 Remote Control | 20 |
| Pro Focus UltraView 2202 System | 20 |
| Battery-Powered Systems | 21 |
| BK Transducers (E14C4t as an example) | 22 |
| Transducer Care and Maintenance | 22 |
| Inspection | 22 |
| Service and Repair | 22 |
| Storing Transducers When Not in Use | 23 |
| Operating and Storage Environment | 24 |
| Covers and Gels | 25 |
| Transducer Covers | 25 |
| Gels | 26 |
| Using Sterile Covers on a System | 26 |
| Chapter 4: Reprocessing Information and Methods | 27 |
| General Information | 28 |
| Levels of Reprocessing | 29 |
| General Precautions | 29 |
| Protecting Transducer Plugs during Immersion | 29 |
| | 30 |
| Reprocessing Equipment | 30 |
| Pre-Cleaning (Point-of-Use Processing) | 31 |
| Leakage Testing | 31 |
| Thorough Cleaning | 32 |
| Manual Disinfection | 33 |
| Reprocessing Other Accessories | 34 |
| Accessories Used for Puncture | 34 |
| Accessories Not Used for Puncture | 34 |
| Cleaning and Disinfecting the System | 35 |
| Reprocessing Remote Controls | 36 |
| Cleaning Remote Controls | 37 |
| Disinfection/Sterilization | 37 |
| Cleaning the Barcode Reader | 38 |
| Automatic Reprocessing Methods | 38 |

| | |
|---|-----------|
| Automated Endoscope Reprocessor – Medivators® Advantage Plus | 38 |
| STERIS SYSTEM 1 Models. | 39 |
| STERRAD Systems | 40 |
| STERIS V-PRO Systems | 42 |
| STERIZONE VP4. | 43 |
| Matachana 130LF, Webeco FA90, Webeco FA95. | 44 |
| Autoclaving. | 46 |
| Chapter 5: Checking and Maintaining Ultrasound Equipment | 47 |
| Check of Equipment Between Each Use | 48 |
| Transducers 20R3, 8838, and 2052. | 49 |
| Transducers 20R3, X14L4, 8838, and 2052 | 49 |
| Transducer 8666-RF, 8809, X18L5s, and I12C4f | 49 |
| Checking the Plug and Watertight Protection Devices Before Immersion. | 49 |
| Checking the Transducer for Leakage | 50 |
| Leakage Testing Kits | 53 |
| Yearly Preventive Maintenance and Performance Test | 53 |
| Yearly Check of Type BF Transducers. | 54 |
| Chapter 6: Watertight Protection Devices | 55 |
| Watertight Plug Lids. | 56 |
| Chapter 7: Disposal | 59 |
| Appendix: Reprocessing Information and Tables | 61 |
| Approved Disinfectants. | 62 |
| Transducer Compatibility | 62 |
| Product and Process Manufacturers | 63 |
| Reprocessing Methods for flex Focus/bk3000/bk3500/bk5000/bkSpecto Systems. | 64 |
| Validated Reprocessing Methods and Material Compatible Reprocessing Methods for the T7P2m (9027) Transducer | 65 |
| Material Compatible Reprocessing Methods for bkFusion Hardware | 66 |
| Validated Reprocessing Methods for the bk3000/bk3500/bk5000/bkSpecto Transducer Series. | 67 |
| Material Compatible Reprocessing Methods for the bk3000/bk3500/bk5000/bkSpecto Transducer Series | 68 |
| Validated Reprocessing Methods for the Flex Focus and Pro Focus Transducer Series | 70 |
| Material Compatible Reprocessing Methods for the Flex Focus and Pro Focus Transducer Series. | 71 |
| Material Compatible Reprocessing Methods for Guides, Attachments and bkFusion Hardware. | 73 |
| Reprocessing Methods for Remote Controls UA1237 and UA2361 | 74 |
| Overview of InstruSafe Instrument Protection Trays for BK Devices. | 75 |
| Reprocessing Table for Craniotomy Transducer 8862 and Burr-Hole Transducer 8863. | 77 |
| Reprocessing Information for Neurosurgical Transducers 8862 & 8863. | 78 |
| Index | 81 |



Chapter 1: Warnings and Cautions









Warnings






| | |
|--|--|
|  Keep plug dry | <p>WARNING T-w5</p> <p>To prevent electrical shock and damage to the transducer, the connector pins in the transducer plug must always be completely dry before you connect to a system.</p> |
|  Transducer, holder and container must be clean | <p>WARNING Reproc-w1</p> <p>To avoid contamination, clean transducers before inserting them into storage containers for transportation.</p> <p>To avoid cross-contamination, make sure that transducer holders and storage containers are clean before inserting clean transducers.</p> |
|  Infection control – follow established procedures | <p>WARNING Reproc-w2</p> <p>Users of this equipment have an obligation and responsibility to provide the highest possible degree of infection control to patients, co-workers and themselves. The instructions in this book are meant as a guide. To avoid cross-contamination, follow all infection control policies (including for reprocessing, packing and storage) for personnel and equipment that have been established for your office, department or hospital.</p> |
|  Pre-clean immediately after use | <p>WARNING Reproc-w3</p> <p>Immediately after use, you must pre-clean the device until visually clean (including any lumens). Conduct the thorough cleaning process as soon as possible after use in order to prevent bioburden drying on the surface. Dried bioburden can lead to inefficient cleaning, disinfection and sterilization, causing a risk of cross-contamination.</p> |
|  Reprocessing | <p>WARNING Reproc-w4</p> <p>Extra care must be taken when cleaning this transducer, as there is no automatic cleaning system that can reprocess the entire transducer.</p> |
|  Screen contamination | <p>WARNING Reproc-w5</p> <p>To avoid contamination of the speaker area, do not touch this area when you turn the monitor to the horizontal or vertical position.</p> |
|  Only approved methods | <p>WARNING Reproc-w6</p> <p>All reprocessing methods described in this book may not apply to all transducers.</p> <p>To prevent damage to a transducer, use only reprocessing methods that are recommended for that transducer.</p> |
|  Neurosurgical covers | <p>WARNING TC-w2</p> <p>For neurosurgical applications, use only non-pyrogenic, sterile probe sheaths (transducer covers) that are approved for neurosurgical use. This means that in the USA they must be market cleared by the FDA and in Europe they must be CE-marked. In Canada, they must be licensed by Health Canada.</p> |
|  | <p>WARNING TC-w7</p> <p>Use only approved non-pyrogenic, sterile transducer covers. This means that in the USA they must be market cleared by the FDA and in Europe they must be CE-marked. In Canada, they must be licensed by Health Canada.</p> |

| | |
|---|---|
|  Creutzfeldt-Jakob | WARNING C-J-w1 Do not use a transducer for neurosurgical applications if the patient is suspected of having Creutzfeldt-Jakob disease. If a neurosurgical transducer has been used on a patient suspected of or diagnosed as being Creutzfeldt-Jakob positive, the transducer must be destroyed, following approved procedures for your hospital. |
|  Damage and reprocessing | WARNING Check-w2 Equipment may be damaged by use or reprocessing. It is important to check it at least once a month (or more often, if it undergoes sterilization) to ensure that it can be effectively reprocessed. If there are any pits or cracks on any equipment surfaces, reprocessing may not give a sterile or disinfected product. |
|  Remote Control | WARNING RC-w1 The remote control requires surface disinfection or sterilization as a minimum. The inside battery compartment cannot be classified as disinfected or sterile. Follow procedures established for your hospital, clinic or institution to avoid cross-contamination when inserting or removing batteries. |
|  Do not use damaged equipment | WARNING Check-w1 To ensure safe operation, do not use the equipment if you find any signs of damage. Contact your BK service representative. If a transducer is dropped, and even if it shows no visible signs of damage, BK recommends that a High Voltage test is conducted before the transducer is used again. |
|  Check of Type BF transducers | WARNING Check-w3 To prevent electrical shock, all transducers that comply with Safety Standard EN60601-1 (IEC60601-1) Type BF must be checked once a year to ensure that they still comply with the requirements of this standard. Transducers that need to be checked have the letters BF or the  symbol printed on them. This check must be carried out by qualified personnel. Contact your BK service representative if you need any help checking your transducers. |
|  Contaminated items | WARNING D-w1 For disposal of contaminated items such as transducer covers or needle guides or other disposable items, follow disposal control policies established for your office, department or hospital. |

Cautions

| | |
|---|--|
|  Physician required in USA | Caution Rx-c1 United States law restricts this device to sale by or on the order of a physician. |
|  Heat, dust, sunlight, condensation | Caution S-c1 To prevent damage to the system, avoid excessive heat, dust and direct sunlight. Also do not use the system if there is visible condensation on it. Wait until it reaches room temperature. |

| | |
|--|--|
|  Handle with care | Caution T-c1 To prevent damage, handle equipment carefully. <ul style="list-style-type: none"> • Don't strike or drop the transducer on a hard surface. Equipment dropped on a hard surface may not be repairable. • Don't step on a cord or run over it with the wheels of the system. |
|  Keep plug dry | Caution Plug-c1 To prevent damage to a transducer or system, protect the plug from contact with liquid. |
|  Immersion: Cover plug – Lid ON | Caution Plug-c2 To prevent damage to the transducer, cover the plug with the watertight protection device before you immerse the transducer and plug in liquid. |
|  Keyboard panel not watertight | Caution S-c2 The keyboard panel of the ultrasound system is not watertight. Be careful not to spill any liquids, gels or moist substances on the keyboard panel. |
|  Gas Processing: Do not cover plug - Lid OFF | Caution Plug-c3 Do NOT use a watertight protection device with any form of gas processing. The transducer can be seriously damaged if a watertight protection device is used. |
|  Do not autoclave transducers | Caution Reproc-c2 Never autoclave the transducers or remote control; this will damage them. |
|  Examine plug and waterproof protection for damage | Caution Plug-c4 Before you reprocess the transducer, inspect the watertight protection device and the transducer plug. If you find any signs of damage, do not immerse the plug. If liquid comes into contact with the plug connector pins, the transducer may be destroyed. Examine the edges of the plug case that contact the lid and also the watertight protection device for cracks and marks. Examine the rubber seal of the plug lid. Look for deep scratches and grooves, holes or tears, brittleness, and looseness anywhere. The transducer or watertight protection device must be checked by a BK service representative if you find signs of damage. |
|  Test for leaks before immersing | Caution Test-c1 You should use the leakage tester to test for leaks. If a transducer is not completely watertight, immersing it can seriously damage it. |












| | |
|--|---|
|  Keep watertight plug lid dry | <p>Caution Test-c4</p> Do not let the watertight plug lid get wet during the testing procedure. Keep it out of the tank. If water gets inside the watertight plug lid, moisture can be transferred from the lid to the plug connector pins during reprocessing. This can damage the transducer. |
|  Do not immerse if pressure drops | <p>Caution Test-c3</p> If the pressure drops to zero after you use the pump, do not place the transducer in the tank. |
|  If you see bubbles, do not release pressure | <p>Caution Test-c5</p> If you see any bubbles, remove the transducer from the tank before you release the pressure. |
|  Do not use test lid for reprocessing | <p>Caution Test-c2</p> The lid part of the leakage tester is for testing only. Do not use it when you reprocess the transducer. |
|  | <p>Caution: T-c5</p> Using a non-recommended disinfection solution, an incorrect solution strength or immersing a transducer longer than recommended by the disinfectant manufacturer can damage the transducer. |













Chapter 2: Markings on the Systems, Transducers and Accessories













Safety Symbols and Information on the Equipment

Table 1 contains brief explanations of the symbols and information used to label the equipment. (Some labels in the table may appear on the transducer.)

The manufacturer disclaims all responsibility for the operating safety, reliability, and performance of the equipment if these symbols and warnings are disregarded in any way.

| Symbol | Name | Description |
|---|---|--|
|  | Medical Device | Device used for medical purpose. |
|  | BK Medical logo | |
|  | Manufacturer label | |
|  | Caution or Warning | Consult accompanying user guide when you encounter this sign on the instrument, to avoid reducing its safety. |
|  | Consult instructions for use | Consult user guide or other instructions. |
|  | CE Marking | Complies with EEC Council Directive 93/42/EEC of 14 June 1993 concerning Medical Devices. |
|  | Rx only | Federal (U.S.A) law restricts sale of this device to physicians or other qualified medical professionals. |
|  | China ROHS 25 Years Lifetime | Environmentally Friendly Use Period for ROHS is 25 years. |
|  | Recycle symbol for corrugated cardboard | Corrugated cardboard for recycling. |
|  | WEEE waste | Within the EU, when you discard waste of electrical and electronic equipment, you must send it to appropriate facilities for recovery and recycling. |
|  | Manufacturer | Legal manufacturer. |

| | | |
|---|-------------------------------|--|
|  | Manufacturer's catalog number | For BK Medical, this is the "Type number" of a product. |
|  | Serial number | Manufacturer's serial number for the specific device. |
|  | Batch code | Manufacturer's batch or lot number for a product. |
|  | Type BF | BF: Isolated from ground Maximum patient leakage current under <ul style="list-style-type: none"> • Normal condition $\leq 100\mu\text{A}$ • Single-fault condition $\leq 500\mu\text{A}$ |
|  | Type BF | BF, defibrillator-proof. |
|  | Type B | B: Maximum patient leakage current under <ul style="list-style-type: none"> • Normal condition $\leq 100\mu\text{A}$ • Single-fault condition $\leq 500\mu\text{A}$ |
| IP57 | Sealing | Dust-protected. Protected against immersion up to 1 m. |
|  | Handle with care | The tip of the transducer is very delicate. Be very careful not to bump the tip. |
|  | Do not reuse | Single-use device. Do not try to process for reuse. Reuse can result in cross-contamination or can compromise the function of the product. |
|  | Do not push | Device may overbalance and fall. |
|  | STERILE | Device is in a sterile condition. |
|  | STERILE EO | Device has been sterilized using ethylene oxide. |
|  | Do not resterilize | Do not resterilize. Resterilization can cause degradation of the materials and compromise the function of the product. |

| | | |
|---|---|--|
|  | Non-sterile | Device is not in a sterile condition. |
|  | Contains latex. | Contains natural rubber latex or latex is present. |
|  | Not made with natural rubber latex | Not made with natural rubber latex. |
|  | Do not use if package or label is damaged | Do not use if product sterilization barrier or its packaging is compromised. |
|  | Fragile, handle with care | Indicates a medical device that can be broken or damaged if not handled carefully. |
|  | Keep dry | Indicates a medical device that needs to be protected from moisture. |
|  | Keep away from sunlight | Indicates a medical device that needs to be protected from light sources. |
|  | Temperature limitation | Keep temperature between the upper and lower limits listed (-20 °C to +60 °C) for BK Medical 88xx/90xx transducers and systems. The 9027 (T7P2m) transducer, BK sterile needle guides, and 3rd party products must be handled according to specific labeling. |
|  | Atmospheric pressure limitation | Indicates the range of atmospheric pressure to which the medical device can be safely exposed. |
|  | Humidity limitation | Keep relative humidity between the upper and lower limits listed. |
|  | Date of manufacture | Date device was manufactured (4 digits for year, 2 digits for month). |
|  | Use by | Last date on which a marked item can be used: expiration date (4 digits for year, 2 digits for month). |








| | | |
|---|---|---|
|  | STERRAD – lid off | Watertight plug lid must not be attached during STERRAD processing. |
|  | Not watertight | Plug must not be immersed. |
|  | Immersion with cap on | Can be immersed if cap is tightened as indicated. |
|  | No immersion with cap off | Must not be immersed if cap is off or not tightened. |
|  | Lock-Unlock | On the watertight plug lid, indicating the locked and unlocked position for the locking pins that hold the lid onto the plug. |
|  | Immersion reprocessing – lid must be on | Watertight plug lid must be attached during immersion. |
|  | Gas reprocessing – lid must be off | Watertight plug lid must not be attached during gas reprocessing or whenever pressure is below 700 hPa (10.2 psi). |

Table 1. Markings on the system, transducers and accessories.

Chapter 3: Introduction & Safety

Introduction


This user guide contains information about caring for and reprocessing BK Medical equipment. It includes important information about what you must do to ensure the safe and proper performance of the system, transducers and accessories. This includes information about cleaning, disinfection and sterilization.


Follow
established
procedures

NOTE: *The instructions below are meant as a guide. They describe the highest level protocol for this level of reprocessing. Always follow the procedures that have been established for your hospital, clinic or institution, as well as any national guidelines.*

Warnings, Cautions, Notes

Pay attention to the difference between Warnings, Cautions and Notes.


| | |
|---|--|
|  | WARNING Warnings contain information that is important for avoiding personal injury. |
|---|--|

| | |
|---|---|
|  | Caution Cautions contain information that is important for avoiding damage to equipment, data |
|---|---|

NOTE: *Notes contain other information that you should be aware of.*

Please find a complete list of warnings and cautions starting on page 6.

General Safety

| | |
|--|---|
|  Physician required in USA | Caution Rx-c1 Federal law in North America restricts this equipment to sale or use by or on the order of a physician. |
|--|---|


bkSpecto



The bk3000/bk5000/bk3500 System and UA2361 Remote Control



The Flex Focus 1202 System and UA1237 Remote Control



Pro Focus UltraView 2202 System



Battery-Powered Systems

The BK battery-powered systems are equipped with high capacity lithium batteries that can power the unit for several hours.

General Recommendations for Battery-Powered Systems

BK recommends that health care professionals and health care facilities take the following steps to help reduce the potential for injury to patients, staff and visitors:


- Do not block any vents.
- Do not tape or attach any object or material to the battery compartment.
- Keep flammable and explosive objects away from battery-powered systems.
- When systems are not in use, but charging, make certain these systems are located in easily visible, fire retardant, non-confined locations away from patient care areas and open sources of oxygen.
- Do not use batteries that do not charge properly. Ensure that batteries are replaced at the recommended replacement intervals. A caution will be displayed on the screen.
- Preventative maintenance information can be found in Chapter 5.
- Notify your BK service representative of damaged batteries.

The FDA has received medical device reports of health hazards associated with batteries used in mobile medical carts and their chargers. These events, which range from smoke production and overheating to equipment fires, can occur with lithium batteries. Note that lithium battery fires are very difficult to extinguish.



In such cases,

- Immediately report the fire according to your hospital protocol. Follow hospital protocol for addressing an energized electrical fire.
- Do not touch the battery.
- Unplug or power off the system if it is safe to do so.
- Remove the system from patient and visitor areas, as safely as possible.

BK Transducers (E14C4t as an example)



Transducer Care and Maintenance

| | |
|---|--|
|  Handle with care | <p>Caution T-c1</p> <p>To prevent damage, handle equipment carefully.</p> <ul style="list-style-type: none"> • Don't strike or drop the transducer on a hard surface. Equipment dropped on a hard surface may not be repairable. • Don't step on a cord or run over it with the wheels of the system. |
|  Keep plug dry | <p>WARNING T-w5</p> <p>To prevent electrical shock and damage to the transducer, the connector pins in the transducer plug must always be completely dry before you connect to a system.</p> |

Inspection

All transducers, the remote control and attachments must be checked regularly to maintain a high level of safety. Details about what to check and when are in the section “Chapter 5: Checking and Maintaining Ultrasound Equipment”.

Service and Repair

If for any reason you must return a transducer to BK Medical, please clean as for storage (if possible). If the transducer is damaged and cannot be reprocessed, clean using disinfection wipes (following the manufacturer’s guidelines) prior to packing and shipping.


Storing Transducers When Not in Use

When storing or transporting a transducer, ensure that:

- The transducer does not get damaged
- The transducer's reprocessing level is maintained

This can be achieved by using a tray with fittings to hold the transducer and the system connector in place (see “Overview of InstruSafe Instrument Protection Trays for BK Devices” on page 75).

If using a container with foam (such as the box the transducer was supplied in), avoid contaminating the foam by ensuring that the transducer is clean, disinfected and dry before placing it in the case. Always reprocess the transducer prior to use.

| | |
|---|--|
|  <p>Transducer, holder and container must be clean</p> | <p>WARNING Reproc-w1</p> <p>To avoid contamination, clean transducers before inserting them into storage containers for transportation.</p> <p>To avoid cross-contamination, make sure that transducer holders and storage containers are clean before inserting clean transducers.</p> |
|---|--|

Storage

- All transducers must at least be thoroughly cleaned and thoroughly dried prior to storage
- The transducers must be stored in a way that maintains reprocessing level (e.g. sterile, high-level disinfected)
- Transducers must be stored in a clean and dry place
- Keep the transducer protected at all times from sharp objects that may damage the transducer or packaging
- The watertight connector lid should not be left on for an extended storage period

See storage and humidity limits for devices on page 25.

Transportation

- In order to prevent damage to the transducer head, lens or connector, use a rigid tray with a lid and internal fittings that keep the transducer and the system connector in place
- Place the transducer in the tray carefully to prevent kinking of the cable
- Before closing the lid, make sure that no part of the transducer is protruding from the tray
- Wrap the tray in plastic material containing air pockets (such as Bubble Wrap), and pack the wrapped tray in a cardboard carton

If using a container with foam (such as the box the transducer was supplied in), avoid contaminating the foam by ensuring that the transducer is clean, disinfected and dry before placing it in the case.

Transducer Holders

Transducer holders can be used for short-term storage of non-critical transducer types, e.g. linear and curvilinear probes. The transducer holders on the system should be cleaned regularly following local guidelines. Semi-critical and critical transducer types, e.g. endocavity and surgical, should be stored in appropriate containers to maintain their level of reprocessing.


Operating and Storage Environment

Systems

Table 2 shows the environmental limits for BK Medical systems during operation and storage.

| | Maximum | Minimum |
|-----------------------|---------------------|--------------------|
| Storage temperature | +60°C (+140°F) | -20°C (-4°F) |
| Operating temperature | +40°C (+104°F) | +10°C (+50°F) |
| Atmospheric pressure | 1060 hPa (15.4 psi) | 700 hPa (10.2 psi) |
| Humidity | 85% RH | |

Table 2. Environmental limits for systems.

| | |
|--|--|
|  Heat, dust, sunlight, conden- sation | Caution S-c1 To prevent damage to the system, avoid excessive heat, dust and direct sunlight. Also do not use the system if there is visible condensation on it. Wait until it reaches room temperature. |
|--|--|

Transducers

Table 3 shows the environmental limits for transducers (and remote controls) during operation and storage.

| | Maximum | Minimum |
|--|-----------------------------|---------------------|
| Storage temperature for: T7P2m | +55°C (+131°F) | -25°C (-13°F) |
| Storage temperature for: all other transducers | +70°C (+158°F) | -25°C (-13°F) |
| Storage humidity | 90% RH | |
| Operating temperature | +40°C (+104°F) | +10°C (+50°F) |
| Operating pressure | 1060 hPa (15.4 psi) | 700 hPa (10.2 psi) |
| Temperature during reprocessing for: T7P2m | +55°C (+131°F) ^a | NA (not applicable) |
| Temperature during reprocessing for: all other transducers | +60°C (+140°F) ^a | NA (not applicable) |
| Pressure during gas processing for: 2052, 8838, 20R3 | NA | 500 hPa (7.3 psi) |
| Pressure during gas processing | NA | 100 hPa (1.5 psi) |

a. Max rate of temperature increase: 15°C/min (27°F/min)

Table 3. Environmental limits for transducers.

Accessories

Follow guidelines for proper storage and handling, as shown by symbols on package label.

Covers and Gels

Transducer Covers

- Transducer covers include sterile and non-sterile condoms and sterile intraoperative transducer and cable covers.
- To reduce the risk of cross-contamination, use a transducer cover when you image.
- You must use a transducer cover for rectal or vaginal imaging. In Germany, you must use a sterile cover for vaginal imaging.
- Follow the procedures that have been established for your hospital, clinic or institution for covers used in conjunction with puncture procedures.

NOTE: *In the United States of America, it is recommended to use probe sheaths (transducer covers) that have been market cleared. In Canada, use only licensed transducer sheaths (covers). In Europe, transducer sheaths must be CE-marked.*

Sterile Covers


BK Medical supplies a range of sterile single-use transducer covers. See the transducer Product Data sheet for appropriate covers for your transducer.


Follow these precautions when putting sterile covers on a transducer:

- Wear sterile gloves.
- When using a puncture attachment, place it gently over the cover and secure it, following the instructions for the puncture attachment (found in the transducer user guide).
- Verify that the cover has not been damaged in the process. If it has, repeat the procedure with a new transducer cover.

Neurosurgical Applications

Special considerations apply to transducer covers for neurosurgical applications.

| | |
|---|---|
|  Neurosurgical covers | WARNING TC-w2 For neurosurgical applications, use only non-pyrogenic, sterile probe sheaths (transducer covers) that are approved for neurosurgical use. This means that in the USA they must be market cleared by the FDA and in Europe they must be CE-marked. In Canada, they must be licensed by Health Canada. |
|---|---|

| | |
|---|--|
|  | WARNING TC-w7 Use only approved non-pyrogenic, sterile transducer covers. This means that in the USA they must be market cleared by the FDA and in Europe they must be CE-marked. In Canada, they must be licensed by Health Canada. |
|---|--|

Gels


Using Sterile Covers on a System


It is important that any cover you use on the system does not affect the readability of the monitor and does not interfere with the touch functionality of the screen or keyboard. Test covers *before* using them during surgical procedures.

Chapter 4: Reprocessing Information and Methods

Personnel must be trained

Personnel in Sterile Processing Departments, as well as nurses, physicians, sonographers and others, may be responsible for reprocessing medical devices. Anyone who reprocesses medical devices should be thoroughly trained in the proper procedures.

| | |
|--|---|
|  Infection control – follow established procedures | WARNING Reproc-w2 Users of this equipment have an obligation and responsibility to provide the highest possible degree of infection control to patients, co-workers and themselves. The instructions in this book are meant as a guide. To avoid cross-contamination, follow all infection control policies (including for reprocessing, packing and storage) for personnel and equipment that have been established for your office, department or hospital. |
|--|---|

| | |
|--|---|
|  Creutzfeldt-Jakob | WARNING C-J-w1 Do not use a transducer for neurosurgical applications if the patient is suspected of having Creutzfeldt-Jakob disease. If a neurosurgical transducer has been used on a patient suspected of or diagnosed as being Creutzfeldt-Jakob positive, the transducer must be destroyed, following approved procedures for your hospital. |
|--|---|


General Information


NOTE: *Reprocessing methods are harsh and can shorten the life of the equipment. Equipment that undergoes reprocessing should be checked regularly. See “Chapter 5: Checking and Maintaining Ultrasound Equipment” on page 47.*

Use the same chemistry each time

BK Medical recommends using the same reprocessing chemistry each time to minimize material degradation caused by chemical interactions from various cleaning, disinfection and sterilization processes.

Reprocessing may cause cosmetic changes to the device material that do not necessarily impact the functionality of the device. The most usual cosmetic changes are color changes in the device materials.

| | |
|--|--|
|  Keep plug dry | WARNING T-w5 To prevent electrical shock and damage to the transducer, the connector pins in the transducer plug must always be completely dry before you connect to a system. |
|--|--|

| | |
|--|---|
|  Damage and reprocessing | WARNING Check-w2 Equipment may be damaged by use or reprocessing. It is important to check it at least once a month (or more often, if they undergo sterilization) to ensure that it can be effectively reprocessed. If there are any pits or cracks on any equipment surfaces, reprocessing may not give a sterile or disinfected product. |
|--|---|

Levels of Reprocessing

The level of processing required depends on the type of equipment and its use.


The CDC (Centers for Disease Control and Prevention) in the USA and the RKI (Robert Koch Institute) in Germany classify medical devices according to their use. For each classification, they specify the level of disinfection/sterilization processing that is required before use. Table 4 summarizes this information.

| Device Classification | Use | Level of Processing Required |
|-----------------------|---|--|
| Noncritical | Device contacts intact skin | Cleaning and disinfection |
| Semi-critical | Device contacts mucous membranes (for example, endocavity applications) | Immediate cleaning and disinfection (in the USA, high-level disinfection or sterilization) |
| Critical | Device enters otherwise sterile tissue (for example, intraoperative applications) Device contacts otherwise sterile tissue (for example, to take a biopsy) | Immediate cleaning followed by sterilization |

Table 4. Levels of disinfection based on device use.

General Precautions


For all types of reprocessing, be sure to observe the environmental limits in Table 3 on page 25. If these limits are exceeded, the transducer may be damaged:

| | |
|--|---|
|  Only approved methods | <p>WARNING Reproc-w6</p> <p>All reprocessing methods described in this book may not apply to all transducers. To prevent damage to a transducer, use only reprocessing methods that are recommended for that transducer.</p> |
|--|---|


Device-specific methods are listed in the reprocessing tables starting on page 64.

If all parts of the complete device are not reprocessed with the same reprocessing steps at the same time, care should be taken not to cross-contaminate the parts of the device, which have already been reprocessed. Some BK transducers are not fully immersible, as the connector cannot be immersed. These parts are therefore to be reprocessed with a wiping method only, and care must be taken to avoid cross-contamination.

Protecting Transducer Plugs during Immersion


| | |
|--|--|
|  Keep plug dry | <p>Caution Plug-c1</p> <p>To prevent damage to a transducer or system, protect the plug from contact with liquid.</p> |
|--|--|

When a transducer is fully immersed (including its plug) during disinfection, the internal components of the plug must *not* get wet. The transducer must be made watertight.

| | |
|---|---|
|  Immersion: Cover plug – Lid ON | Caution Plug-c2 To prevent damage to the transducer, cover the plug with the watertight protection device before you immerse the transducer and plug in liquid. |
|---|---|

Reprocessing Equipment

Proper cleaning is essential to the success of any disinfection or sterilization procedure. Equipment must be cleaned immediately after it is used and before it is disinfected or sterilized.

| | |
|--|---|
|  Pre-clean immediately after use | WARNING Reproc-w3 Immediately after use, you must pre-clean the device until visually clean (including any lumens). Conduct the thorough cleaning process as soon as possible after use in order to prevent bioburden drying on the surface. Dried bioburden can lead to inefficient cleaning, disinfection and sterilization, causing a risk of cross-contamination. |
|--|---|

Before cleaning, always remove covers, accessories and attachments, including dummy attachments. Then clean the transducer and reusable attachments thoroughly. See the cleaning steps listed in the following pages.




Figure 1. Overview of reprocessing steps.

The general cleaning process should include the following overall steps, performed as efficiently as possible. You must also follow your local and/or national guidelines.

- 1 Pre-cleaning at point of use
- 2 Leak testing

- 3 Thorough cleaning including rinsing or processing through an automated washer system
- 4 Disinfection or sterilization

Pre-Cleaning (Point-of-Use Processing)

For pre-cleaning of equipment, follow your local and/or national guidelines.

NOTE: *If you are pre-cleaning a transducer and immersing or rinsing it in step 2, 3 and/or 4, test the transducer for leakage after you unplug it and before you put it in liquid (immersion or rinsing with running water). The leakage testing procedure is described on page 52.*

Procedure for pre-cleaning:

- 1 Unplug the transducer from the system.
- 2 Immediately remove any cover, puncture guides or other attachments and disassemble all parts. Wipe off any gel or biological material with a moist cloth or sponge moistened with detergent solution (made using manufacturer's instructions) or water, taking care to remove all visible contamination.
Note that the cloth or sponge should be disposed of, sterilized, or high-level disinfected after each use.
- 3 If the device has a lumen, make sure that you brush it with a moistened brush (brush size compatible with the lumen) until the lumen is visually clean.
- 4 Do not attempt to put a brush or anything else through a water inlet. If you suspect that a water channel has become contaminated, flush it with detergent.
- 5 If you have used a detergent solution, follow the manufacturer's instructions regarding rinsing and wiping.
- 6 If transportation is required, transport the pre-cleaned transducer to the reprocessing area in a closed container that prevents exposure of staff, patients, or the environment to potentially infectious materials. Containers, sinks, and basins should be large enough that the transducer will not be damaged.

Leakage Testing

For a full description of the leakage tester, see "Leakage Testing Procedure" on page 52.

Thorough Cleaning

Validated detergents are listed in the tables starting on page 73. Dilute and use according to the detergent manufacturer's instructions.

Manual Cleaning by Immersion

Procedure for manual cleaning by immersion:

- 1 Fill a sink or bowl with freshly-made detergent solution.
- 2 Immerse the transducer and all removable parts, and, if possible, the cable and plug¹. Make sure that the watertight plug lid is attached correctly and the transducer has been tested for leakage. See “Leakage Testing Procedure” on page 52.

NOTE: *To prevent the splashing of contaminated fluid and aerosolization of bioburden, the device should be submerged in the detergent solution when you perform all subsequent cleaning steps.*

- 3 Use a suitable soft brush (for example, a soft nail sponge/brush like surgeons use) to thoroughly clean all parts of the device, paying special attention to the tip, any lumens, buttons, lever, edges or grooves.
- 4 Do not attempt to put a brush or anything else through a water inlet. If you suspect that a water channel has become contaminated, flush it with the detergent.
- 5 To clean device lumens (e.g. in transducers or biopsy needle guides), use a brush appropriate to the size of the lumen. Brush the lumen thoroughly with a circular movement for about 10 seconds (at least three times); make sure that all inner surfaces of the lumen have been brushed. Repeat from the other end of the lumen. After each passage, remove any visible debris from the brush before reinserting it.

NOTE: *When cleaning the transducer lumen, rub gently with the brush to avoid causing damage.*

- 6 Flush all channels with the detergent solution to remove debris. If any debris is noted in the flushing solution repeat the previous step.
- 7 Continue soaking the transducer, puncture attachments, and any internal lumens until the total detergent contact time specified by the manufacturer is reached.
- 8 Visually inspect for any remaining soil and if necessary repeat the steps, starting at step 3.
- 9 Make sure to follow the rinse (and neutralization) instructions from the detergent manufacturer for rinsing the device and all removable parts.
If insufficient instructions are supplied, then, as a minimum, thoroughly rinse the device and all removable parts with running water (preferably deionized or distilled or RO water (purified by reverse osmosis)) with a flow rate of approximately 2 l/min and a temperature between 10°C (50°F) and 40°C (104°F) until all signs of residual debris and cleaning solution are removed (for approximately 1 min).

1. Some transducer types are not provided with a plug lid, and therefore the connector must not be immersed. See “Transducers Excluded from Leakage Testing” on page 51.

Ensure that any lumens, buttons, lever, edges or grooves are thoroughly rinsed and that the lumen is flushed with water at least 2 times.

- 10** Remove water from all lumens and the exterior of the device with a clean disposable soft, lint-free cloth or air dry or use a drying cabinet (do not exceed 60°C (140°F)).

The device is now ready for disinfection or sterilization.

Manual Cleaning by Wiping

For manual cleaning by wiping, follow wipe manufacturer's instructions. Make sure to clean all surfaces and comply with the prescribed contact time.

Manual Disinfection

Before you disinfect a transducer, it must be thoroughly cleaned. This includes being pre-cleaned at point of use and then cleaned, rinsed and dried. See "Reprocessing Equipment", starting on page 30, for more information. Validated and material compatible reprocessing methods are listed in the tables starting on page 73.

Manual Disinfection by Immersion

To disinfect by immersion:

- 1** Use a disinfectant method that has been approved (or evaluated for material compatibility) for the transducer. See "Appendix: Reprocessing Information and Tables" on page 61.
- 2** Follow the disinfectant manufacturer's instructions for procedure and immersion times.

Make sure that the solution passes through any built-in lumens or grooves. If necessary, use a suitable brush to make sure there are no air bubbles in the channel.

- 3** If specified by the disinfectant manufacturer, rinse off the disinfectant thoroughly with water (preferably deionized or distilled or RO (purified by reverse osmosis) water) with a temperature between 10°C (50°F) and 40°C (104°F), thoroughly flushing any channels. Follow the disinfectant manufacturer's instructions for procedure and volume of water.

For high-level disinfectants, the water should be sterile, and you should wear sterile gloves.

- 4** Dry with a clean soft cloth (sterile for high-level disinfection).
- 5** Examine the transducer and the cable and connector for signs of damage.

Manual Disinfection by Wiping

For manual disinfection by wiping, follow wipe manufacturer's instructions. Make sure to disinfect all surfaces and comply with the prescribed contact time.

Reprocessing Other Accessories

Unless this guide contains alternative instructions, follow the manufacturer's instructions for cleaning and disinfecting any accessories such as movers and steppers and brachy matrices.

Accessories Used for Puncture

All puncture attachments must be pre-cleaned, cleaned and sterilized before use, unless they are supplied sterile. See "Reprocessing Equipment", starting on page 30.

Accessories Not Used for Puncture

Accessories that are not used for puncture are classified as semi-critical or non-critical devices. That means that they may not need to be sterilized or autoclaved and that cleaning and disinfection may be enough. See "Reprocessing Equipment", starting on page 30 for more information.

NOTE: *Accessories that contact mucous membranes (for example, devices used in body cavities) are classified as semi-critical, and in the USA, the reprocessing level must be high-level disinfection or sterilization. For disinfection by immersion, follow the procedures in "Manual Disinfection by Immersion" on page 33.*

Dummy Channel Brackets UA1272 and UA1325-w

The dummy channel brackets can be cleaned in the same way as described for attachments. See "Reprocessing Equipment", starting on page 30

After pre-cleaning and cleaning, the brackets can be autoclaved, as described on page 46.

Reusable Water Standoff System for Anorectal 3D Transducer and 20R3 Transducer

The following parts of the water standoff system must be cleaned as described in "Reprocessing Equipment" on page 30:

- water standoff collar (UA0671)
- rectosigmoidoscope and obturator (UA0672 or UA0673)
- O-rings (UA0674)

The water standoff collar, O-rings, rectosigmoidoscope and obturator can be autoclaved.

NOTE: *Separate the O-rings from the water standoff collar before cleaning and disinfecting as well as before autoclaving.*

Magnetic Wheel Mover UA 0513


Please refer to the magnetic wheel mover user guide for cleaning and disinfection instructions.

Cleaning and Disinfecting the System

Precautions

Although the system surface is resistant to chemicals, strong chemicals may discolor it.

The control panel is sealed underneath and is designed to resist limited amounts of liquid. Liquids should not be poured on it, however. See the caution below.

| | |
|--|--|
|  <p>Keyboard panel not watertight</p> | <p>Caution S-c2</p> <p>The keyboard panel of the ultrasound system is not watertight. Be careful not to spill any liquids, gels or moist substances on the keyboard panel.</p> |
|--|--|

Always turn off the system power before cleaning. If possible, disconnect or unplug the power cord.

Follow the manufacturer's instructions for any cleaning and disinfection products you use.

Cleaning the System

See a list of validated and material compatible system cleaning products on page 64.

Clean the system, including the hand rest, transducer holders (including endo transducer holders) and keyboard panel (and touch screen, if relevant) after every examination. Do not let biological material dry on the system. Clean the monitor if it has been touched.

To clean the system cabinet (including battery compartment):

- 1 If your system has a hand rest, remove it and reprocess it separately. You can also remove transducer holders and the gel holder for reprocessing.
- 2 Use a soft non-abrasive cloth moistened with a mild, general purpose, non-abrasive detergent solution - or use a wipe product manufactured for this purpose. For more information, see "Reprocessing Methods for Remote Controls UA1237 and UA2361" on page 74, or see the system Product Data sheet.
- 3 Wipe the system.
- 4 If necessary, use a damp cloth to remove any detergent residue.
- 5 Wipe dry with a lint-free cloth.
- 6 The trackball can also be removed for cleaning. See below.

To clean the monitor and/or touch screen:

- 1 Use a soft cloth and, if necessary, a product listed in the system Product Data sheet.
- 2 Gently wipe the monitor face/touch screen.
Make sure not to scratch the monitor/touch screen.

To clean the control panel:

- 1 Moisten a soft, non-abrasive cloth with a mild, general purpose, non-abrasive detergent solution - or use a wipe product manufactured for this purpose.
- 2 Wipe the control panel.
- 3 Use a cotton swab to clean around keys or controls. Use a toothpick to gently remove solids from between keys and controls.

When cleaning the operator control panel, make sure not to spill or spray any liquid on the controls, into the system cabinet, or in the transducer sockets.

To clean the trackball:

The trackball can be removed completely for cleaning and disinfection.

- To remove the trackball, rotate the ring around the trackball counterclockwise and lift it off. Remove the trackball.
- After you replace the trackball, replace the ring and rotate it clockwise to tighten it.


Disinfecting the System

- All parts of the system, including the monitor and battery compartment, can be wiped down with validated disinfectants. For more information, see “Reprocessing Methods for Remote Controls UA1237 and UA2361” on page 74, or see the system Product Data sheet.

Reprocessing Remote Controls

The remote control itself is a sealed unit and can be totally immersed if the battery cap is screwed on tight.


The interior of the remote cannot be cleaned and is therefore not to be considered cleaned, disinfected or sterilized. Take care when inserting batteries.

| | |
|---|---|
|  Remote Control | <p>WARNING RC-w1</p> <p>The remote control requires surface disinfection or sterilization as a minimum. The inside battery compartment cannot be classified as disinfected or sterile. Follow procedures established for your hospital, clinic or institution to avoid cross-contamination when inserting or removing batteries.</p> |
|---|---|

For validated reprocessing methods for the remote controls, see “Reprocessing Methods for Remote Controls UA1237 and UA2361” on page 74.

Before cleaning or immersing in disinfectant (including STERIS SYSTEM 1 processing):

- Screw the battery cap on tight until the arrow points to the area of the battery cap with a large gap between the ridges.




Cleaning Remote Controls

Proper cleaning is essential for the success of any disinfection or sterilization procedure. Equipment must be cleaned immediately after it is used and before it is disinfected or sterilized. Use a method that has been approved (or evaluated for material compatibility) for the remote control. See “Reprocessing Methods for Remote Controls UA1237 and UA2361” on page 74.

To clean remote controls (surface cleaning only):

- 1 Make sure that battery cap is screwed on tight so that the arrow points to the area of the battery cap with a large gap between the ridges.
- 2 Immediately rinse or wipe off any visible contamination (such as biological substances) with a detergent solution or water between 10 °C (50 °F) and 40 °C (104 °F), using a brush if necessary.
- 3 Clean with a detergent (see “Reprocessing Methods for Remote Controls UA1237 and UA2361” on page 74) and a soft-bristled nail brush (like surgeons use) or cloth to remove proteins/soil. Follow detergent manufacturer’s guidelines.
- 4 Rinse thoroughly with running tap water between 10 °C (50 °F) and 40 °C (104 °F).
- 5 Dry with a disposable cloth or air dry.
- 6 Thoroughly examine all surfaces that have been cleaned and visually inspect the entire device to make sure it is clean.



Disinfection/Sterilization

Start by cleaning (following recommend steps above).


After the remote control has been cleaned, it can be disinfected or sterilized. See “Reprocessing Methods for Remote Controls UA1237 and UA2361” on page 74 for more information.

Before and after each processing, you must check remote controls for surface pits and cracks.

Before you put a remote control into a STERRAD or V-Pro System:

- Unscrew the cap and remove the batteries. Leave the cap off.
- Reprocess the batteries *with* the remote control and cap but not *in* it.

NOTE: You must use VARTA brand alkaline batteries (LR6, AA) if you process the remote control in a STERRAD system.



Cleaning the Barcode Reader

To clean the barcode reader:

The barcode reader's design allows safe cleaning of the product with a variety of cleaning products and disinfectants. The manufacturer recommends wiping the barcode reader with the following list of approved cleaners:

- Isopropyl alcohol (suitable for medical devices)
- Bleach/sodium hypochlorite (suitable for medical devices)
- Hydrogen Peroxide
- Gentle dish soap and water

Automatic Reprocessing Methods

Proper cleaning is essential for the success of any disinfection or sterilization procedure. All covers and attachments must be removed and all channels thoroughly cleaned. The equipment must be cleaned immediately after it is used and *before* it is reprocessed further.

NOTE: *The following validated automated reprocessing methods state the number of reprocessing cycles used in the validation. The validated number of cycles is not the maximum number of possible cycles for the device, unless this is specifically stated. In general, the maximum number of cycles cannot be stated, as this depends on the usage and handling of the device (e.g. use, all reprocessing steps, storage and transport). By thoroughly performing the checks described in "Chapter 5: Checking and Maintaining Ultrasound Equipment" on page 47, any signs of damage will determine the maximum number of cycles.*

Automated Endoscope Reprocessor – Medivators® Advantage Plus

High-level disinfection

Medivators® Advantage Plus is a high level disinfection reprocessing system for cleaned immersible, reusable, heat-sensitive medical devices. The reprocessing cycle includes a washing step, which in the USA has FDA clearance for pre-cleaned devices.


A number of BK Medical transducers have been validated in Medivators Advantage Plus with Medivators Intercept (detergent), Medivators Rapiocide PA High-Level Disinfectant and flushing with 70% isopropyl alcohol.

Number of Cycles Validated for Transducers

Medivators AER reprocessing has been validated for 100 cycles.

Recommended Cycles

- Medivators Advantage Plus Endoscope Reprocessing System: Lid on

| | |
|---|--|
|  <p>Immersion: Cover plug – Lid ON</p> | <p>Caution Plug-c2</p> <p>To prevent damage to the transducer, cover the plug with the watertight protection device before you immerse the transducer and plug in liquid.</p> |
|---|--|

Follow the manufacturer's instructions for using Medivators AER systems.

STERIS SYSTEM 1 Models

Liquid
chemical
sterilization

STERIS SYSTEM 1 models (STERIS SYSTEM 1¹, SYSTEM 1E, SYSTEM 1 PLUS¹ and SYSTEM 1 EXPRESS¹) are low temperature systems for liquid chemical sterilization of cleaned, immersible, and reusable critical and semi-critical heat-sensitive medical devices. These processes involve immersing items in a circulated solution containing a sterilant concentrate (S40™) that contains peracetic acid, followed by rinses.

A number of BK Medical transducers have been designed to be compatible with STERIS SYSTEM 1 models. The transducer plugs must be fitted with a watertight protection device during processing with STERIS systems. See “Chapter 6: Watertight Protection Devices” on page 55.

Number of Cycles Validated

STERIS SYSTEM 1 models have been validated for 100 cycles.

Recommended Cycles

The following cycles are recommended for BK Medical devices:

Transducers

| | |
|----------------------------|--------|
| • STERIS SYSTEM 1: | Lid on |
| • STERIS SYSTEM 1E: | Lid on |
| • STERIS SYSTEM 1 Plus: | Lid on |
| • STERIS SYSTEM 1 Express: | Lid on |

Remote Control (surface sterilization)

| | |
|----------------------------|--------------------------|
| • STERIS SYSTEM 1: | Batteries inside, lid on |
| • STERIS SYSTEM 1E: | Batteries inside, lid on |
| • STERIS SYSTEM 1 Plus: | Batteries inside, lid on |
| • STERIS SYSTEM 1 Express: | Batteries inside, lid on |

Follow the manufacturer’s instructions for using STERIS systems. Follow the information from the manufacturer of the system about the correct use of the tray for each version of the system 1 (can be found on www.steris.com).

STERIS Quick Connect Table:

| Transducer | SYSTEM 1 - rest of world | SYSTEM 1E - USA only | SYSTEM 1 Plus & 1 Express |
|------------|---------------------------|-------------------------|---------------------------|
| 8808 | Quick Connect: QMC1733INT | Quick Connect: QMC1733E | Quick Connect: QMC1733INT |
| | Tray: C1220S1 | Tray: C1220E | Tray: C1220INT |

1. STERIS SYSTEM 1, SYSTEM 1 Plus and SYSTEM 1 Express have not been market cleared by the FDA.

STERRAD Systems

Liquid
chemical
sterilization

STERRAD 100S, 200, STERRAD NX, STERRAD 100NX

STERRAD® systems use low-moisture and low-temperature hydrogen peroxide gas plasma technology to process heat- and moisture-sensitive medical devices. The systems are intended for terminal sterilization of properly cleaned, rinsed, and thoroughly dried reusable medical devices.

Number of Cycles Validated for Transducers

STERRAD system processing has been validated for 100 cycles.

Number of Cycles Validated for Remote Control

Because of material degeneration, STERRAD system processing has been validated for a maximum of 50 processing cycles for the remote control.

Recommended Cycles


The following cycles are recommended for BK Medical devices:

Transducers

| | | |
|---------------|----------------|----------------------------------|
| STERRAD 100NX | Standard cycle | Lid off |
| | Express cycle | Lid off |
| STERRAD NX | Standard cycle | Lid off |
| STERRAD 100S | Short cycle | Lid off (for US: one cycle only) |
| STERRAD 200 | Short cycle | Lid off |

Remote Control (surface sterilization). Batteries and lid to be sterilized separately.

| | | |
|----------------------|----------------|----------------------------------|
| STERRAD NX and 100NX | Standard cycle | Lid off |
| STERRAD 100S | Short cycle | Lid off (for US: one cycle only) |
| STERRAD 200 | Short cycle | Lid off |

| | |
|--|--|
|  Gas Processing: Do not cover plug - Lid OFF | <p>Caution Plug-c3</p> <p>Do NOT use a watertight protection device with any form of gas processing. The transducer can be seriously damaged if a watertight protection device is used.</p> |
|--|--|

Follow the manufacturer's instructions for using STERRAD systems, including instructions for packaging devices before processing them.

APTIMAX trays are indicated for use in the STERRAD Sterilization System. Instrusafe trays (see “Overview of InstruSafe Instrument Protection Trays for BK Devices” on page 75) have been approved for STERRAD processing of BK devices.

NOTE: *When using a tray, ensure that no parts of the device are pressed hard against the sides of the tray.*

STERIS V-PRO Systems

V-PRO maX, V-PRO 1 Plus, V-PRO 60 and V-PRO 1

V-PRO systems use low-moisture and low-temperature vaporized hydrogen peroxide technology to process heat- and moisture-sensitive medical devices. They are intended for terminal sterilization of properly cleaned, rinsed, and thoroughly dried reusable medical devices.

A number of BK Medical transducers have been designed to be compatible with the V-PRO systems.

Recommended Cycles

The following cycles are recommended for BK Medical devices:

Transducers

| | | |
|--------------|-----------------------------------|---------|
| V-PRO maX | Non lumen cycle or Flexible cycle | Lid off |
| V-PRO 1 Plus | Non lumen cycle | Lid off |
| V-PRO 60 | Non lumen cycle | Lid off |

Remote Control (batteries and lid to be sterilized separately)

| | | |
|--------------|----------------|---------|
| V-PRO maX | Lumen cycle | Lid off |
| V-PRO 1 Plus | Lumen cycle | Lid off |
| V-PRO 60 | Lumen cycle | Lid off |
| V-PRO 1 | Standard cycle | Lid off |

Number of Cycles Validated: STERIS V-PRO systems reprocessing has been validated for 100 cycles.

| | |
|--|--|
| Gas Processing: Do not cover plug - Lid OFF | <p>Caution Plug-c3</p> <p>Do NOT use a watertight protection device with any form of gas processing. The transducer can be seriously damaged if a watertight protection device is used.</p> |
|--|--|

Follow the manufacturer's instructions for using STERIS V-PRO systems, including instructions for packaging devices before processing them.

STERIS trays no. VP0040, VP0041, VP0042, VP0043 and VP0044 have been approved for V-PRO processing of BK devices.

Instrusafe trays (see "Overview of InstruSafe Instrument Protection Trays for BK Devices" on page 75) have been approved for V-PRO processing of BK devices.

NOTE: *When using a tray, ensure that no parts of the device are pressed hard against the sides of the tray.*

STERIZONE VP4

STERIZONE[®] VP4 sterilizer uses vaporized hydrogen peroxide (H₂O₂) and ozone (O₃) in a multiphase process to sterilize heat-sensitive medical devices. The sterilizer is intended for terminal sterilization of properly cleaned, rinsed, and thoroughly dried reusable medical devices.

A number of BK Medical transducers have been validated to be compatible with the STERIZONE VP4 sterilizer. Please see the list on page 67 and page 68.

Recommended Cycles

The sterilizer offers a single preset sterilization cycle (Cycle 1) designed for the sterilization of a wide variety of loads.

Transducers

| | | |
|---------------|-------------------------|---------|
| STERIZONE VP4 | Only one possible cycle | Lid off |
|---------------|-------------------------|---------|

Number of Cycles Validated: STERIZONE VP4 processing has been validated for 100 cycles.

| | |
|--|---|
| Gas Processing: Do not cover plug - Lid OFF | Caution Plug-c3 Do NOT use a watertight protection device with any form of gas processing. The transducer can be seriously damaged if a watertight protection device is used. |
|--|---|

Follow the manufacturer's instructions for using STERIZONE VP4 systems, including instructions for packaging devices before processing them.

NOTE: *When using a tray, ensure that no parts of the device are pressed hard against the sides of the tray.*

Matachana 130LF, Webeco FA90, Webeco FA95

Low Temperature Steam Formaldehyde Systems (LTSF)

The LTSF systems are formaldehyde sterilizers using 2% formaldehyde for sterilization of medical devices. A number of BK Medical transducers have been designed to be compatible with these systems. Please see the list on page 67 and page 68.

Recommended Cycles



Always use the “60°C sterilization program” cycle for sterilization of approved BK Medical devices:

Transducer Sterilization:

| System: | Cycle: | LTSF Valve: |
|-----------------|-----------------------------------|---|
| Matachana 130LF | Cycle: 60°C sterilization program | EasyENDOVALVE: Matachana art.-no 76330 (or previously Webeco art.-no. 046763) |
| Webeco FA90 | Cycle: 60°C sterilization program | |
| Webeco FA95 | Cycle: 60°C sterilization program | |

Transducer Series: Please see the list on page 67 and page 68.

Connector lid and LTSF valve must be attached during processing.

| | |
|---------------------------|--|
| bk3000/bk3500/bk5000: | <p>BK connector lid UA1418 with Matachana art.-no 76330 valve</p>  |
| Flex Focus and Pro Focus: | <p>BK connector lid UA1419 with Matachana art.-no 76330 valve</p>  |

Number of Cycles Validated: The LTSF systems have been validated for 100 cycles.


Prior to sterilization, attach the LTSF valve to connector lid UA1418 or UA1419. When using a sterile barrier system (e.g. wrapping), the transducer and valve must be placed into the sterile barrier system together.

After sterilization, the transducer and valve must remain wrapped for transport and storage. Immediately before use, unwrap the device, aerate the inside of connector lid by turning the screw on the valve, and unscrew the LTSF valve from the connector. Remove the connector lid and valve and place them in a protective bag.

Follow the manufacturer's instructions for using and reprocessing LTSF systems and the LTSF valve.

Autoclaving

All BK Medical steel puncture attachments can be autoclaved after cleaning. (See “Reprocessing Equipment” on page 30.) Autoclaving sterilizes using steam (from water) under pressure. Other transducer accessories may be autoclavable (contact your local BK representative for information). See “Material Compatible Reprocessing Methods for Guides, Attachments and bkFusion Hardware” on page 73.

| | |
|---|---|
|  Do not autoclave transducers | Caution Reproc-c2 Never autoclave the transducers or remote control; this will damage them. |
|---|---|

To sterilize steel parts by autoclaving:

- | | |
|---------------------------|--|
| Packaging for autoclaving | 1 Pack all cleaned parts in a pouch suitable for steam sterilization, or in a tray with paper wrap according to EN ISO 11607 “Packaging for Terminally Sterilized Devices”, or according to local hospital procedure. Follow the pouch manufacturer’s specifications or the local regulations for how to pack and seal the pouches. |
| Autoclaving, parameters | 2 Autoclave all parts of the puncture attachment or accessory, including reusable needles and needle guides. The suggested process parameters for sterilizing attachments are the following: <ul style="list-style-type: none"> • 4 vacuum cycles 100–1000 hPa (1.5–14.5 psi) • Sterilization cycle 134°C (273°F) for 3 min • Cooling phase 100 hPa (1.5 psi) for 5 min |

To sterilize non-steel parts by autoclaving:

For non-steel parts that can be autoclaved, follow the above steps, but be aware that exceeded temperatures can deform (bend) and otherwise damage the plastic parts. See page 73 for information on which non-steel parts cannot be autoclaved, and how to reprocess parts that cannot be autoclaved.

Chapter 5: Checking and Maintaining Ultrasound Equipment

Ultrasound equipment requires regular checks and maintenance. Table 5 contains a list of required checks.

NOTE: *If you find any cracks or irregularities on the transducer, please contact your local BK Medical representative.*


| What to check | When to check |
|--|---|
| Transducer, connector, cable, remote control, attachments and reusable puncture guides for cracks and irregularities See: "Check of Equipment Between Each Use" on page 48. | Between each use |
| Transducer for leakage See: "Checking the Transducer for Leakage" on page 50. | See recommended frequency on page 50 |
| Watertight protection device and transducer plug's waterproof gaskets and seal for cracks and marks See: "Checking the Plug and Watertight Protection Devices Before Immersion" on page 49. | Before immersing the transducer and/or the plug for cleaning or disinfection OR at least monthly (or more often in cases of heavy use) |
| Preventive maintenance and performance test of entire system See: "Yearly Preventive Maintenance and Performance Test" on page 53. | Yearly |
| Type BF transducers to make sure they still comply with requirements See: "Yearly Check of Type BF Transducers" on page 54. | |

Table 5. Required checks of ultrasound equipment.

Check of Equipment Between Each Use

For reprocessing methods to be effective, external surfaces must be in good condition. Transducers (including connector, cable, remote control, attachments and reusable puncture guides) should be checked between each use for signs of damage.


How often For critical-use and semi-critical applications, you should carry out a detailed check for damage between each use. For non-critical applications, this check can be performed daily. For a list of applications, please see the Recommended Leakage Testing Frequency table on page 51.

| | |
|--|--|
|  Damage and reprocessing | <p>WARNING Check-w2</p> <p>Equipment may be damaged by use or reprocessing. It is important to check it at least once a month (or more often, if it undergoes sterilization) to ensure that it can be effectively reprocessed. If there are any pits or cracks on any equipment surfaces, reprocessing may not give a sterile or disinfected product.</p> |
|--|--|

- Damage signs
- Check the transducer for the following signs of damage:
- Pits or cracks anywhere
 - Deep scratches on any surfaces
 - Splitting or peeling of the sealant around the front face (acoustic surface)
 - Damage to the joint filler on the body of the transducer
 - Damage to, or evidence of, contamination on the pins of the transducer plug
 - Damage to the cable or cable bonding around the cable flex relief
 - Deformation or other damage (to non-steel puncture guides) caused, for example, by autoclaving with excessive heat
 - Blistering on the surface of the remote control

NOTE: *The front face (acoustic surface) must appear uniform and be fully attached to the rest of the transducer. It must not be swollen or peeling off*

NOTE: *Use good light and a magnifier when checking the transducer*

| | |
|---|--|
|  <p>Do not use damaged equipment</p> | <p>WARNING Check-w1</p> <p>To ensure safe operation, do not use the equipment if you find any signs of damage. Contact your BK service representative.</p> <p>If a transducer is dropped, and even if it shows no visible signs of damage, BK recommends that a High Voltage test is conducted before the transducer is used again.</p> |
|---|--|

Transducers 20R3, 8838, and 2052

Before you use one of these transducers, thoroughly inspect the connector plug and the rubber gasket on the plug (this refers to the small, round plug).

Transducers 20R3, X14L4, 8838, and 2052

Damage to the surface of these transducers may lead to oil leaking from the transducer.

Transducer 8666-RF, 8809, X18L5s, and I12C4f


Before you use one of these transducers, inspect the flexible, black rubber next to the acoustic array, and the articulation joint on the flexible tip of the transducer to make sure that there are no defects in these areas.

Checking the Plug and Watertight Protection Devices Before Immersion


To make sure that liquid does not get into a plug during immersion, the watertight protection device must be dry inside and it must make a tight seal with the plug¹.

Inspect the equipment for signs of damage to the plug, lid and the rubber sealing.


1. Some transducer types are not provided with a plug lid, and therefore the connector must not be immersed. See "Transducers Excluded from Leakage Testing" on page 51.

| | |
|--|---|
|  Examine plug and waterproof protection for damage | <p>Caution Plug-c4</p> <p>Before you reprocess the transducer, inspect the watertight protection device and the transducer plug. If you find any signs of damage, do not immerse the plug. If liquid comes into contact with the plug connector pins, the transducer may be destroyed.</p> <p>Examine the edges of the plug case that contact the lid and also the watertight protection device for cracks and marks. Examine the rubber seal of the plug lid. Look for deep scratches and grooves, holes or tears, brittleness, and looseness anywhere.</p> <p>The transducer or watertight protection device must be checked by a BK service representative if you find signs of damage.</p> |
|--|---|


Checking the Transducer for Leakage

| | |
|---|--|
|  Do not use damaged equipment | <p>WARNING Check-w1</p> <p>To ensure safe operation, do not use the equipment if you find any signs of damage. Contact your BK service representative.</p> <p>If a transducer is dropped, and even if it shows no visible signs of damage, BK recommends that a High Voltage test is conducted before the transducer is used again.</p> |
|---|--|

Before you immerse a transducer, check the transducer for cracks and irregularities. See “Check of Equipment Between Each Use” on page 48 for more information.

| | |
|---|--|
|  Immersion: Cover plug – Lid ON | <p>Caution Plug-c2</p> <p>To prevent damage to the transducer, cover the plug with the watertight protection device before you immerse the transducer and plug in liquid.</p> |
|---|--|

If the transducer is fitted with a watertight lid, it is recommended to use the appropriate leakage tester (UA1414 or UA1404) to verify that the transducer is watertight.

| | |
|--|---|
|  Test for leaks before immersing | <p>Caution Test-c1</p> <p>You should use the leakage tester to test for leaks. If a transducer is not completely watertight, immersing it can seriously damage it.</p> |
|--|---|

Recommended Leakage Testing Frequency

Check transducers for leakage regularly, as a leakage may impair the performance and safety of the equipment. Recommended leakage testing frequency depends on the design and use of the transducer.

Transducers Excluded from Leakage Testing

| | |
|-------|---|
| 8819 | Excluded due to transducer design or no immersible plug lid |
| 8830 | |
| 8837 | |
| 8670 | |
| 2052 | |
| 8838 | |
| 5C1e | |
| 14L3e | |
| 5P1e | |

Leakage Testing Table


| Transducer Classification | Applications | Before rinse and/or immersion | Leakage testing after each use | Minimum leakage testing frequency |
|---------------------------|--|-------------------------------|--------------------------------|-----------------------------------|
| Non-critical use | Abdominal Adult Cephalic (Transcranial) Cardiac Adult Fetal, including Obstetrics Musculoskeletal Vessel (Peripheral Vessel) Small Parts (also called Small Organs) Neonatal Cephalic | Always | No | Recommended every 3 months |
| Semi-critical use | Transrectal Transvaginal | Always | No | Recommended every 3 months |
| Critical use | Intraoperative Intraoperative (Neuro) | Always | Yes | Recommended after each use |

Leakage Testing Setup

After the transducer plug is covered with a special test lid, air is pumped into the transducer. The transducer and covered plug are then placed in a tank filled with water. If bubbles appear, it is a sign that the transducer, cable, or plug contains a hole and is not watertight.




Figure 2. Example of a leakage testing setup with UA1414. Look for bubbles in the water.


| | |
|---|--|
|  Keep watertight plug lid dry | <p>Caution Test-c4</p> <p>Do not let the watertight plug lid get wet during the testing procedure. Keep it out of the tank.</p> <p>If water gets inside the watertight plug lid, moisture can be transferred from the lid to the plug connector pins during reprocessing. This can damage the transducer.</p> |
|---|--|

Leakage Testing Procedure

To test a transducer for leaks:


- 1** Visually check the edge of the plug and the seal on the lid to make sure that they will fit tightly. Place the test lid on top of the plug with the locking pins unlocked (pointing at the open lock sign). Make sure that the lid is properly aligned and seated on top of the plug.
- 2** Firmly attach the test lid, as described in “Watertight Plug Lids” on page 56.
- 3** Pump slowly to increase the pressure up to a relatively stable level of 150 mm Hg. Keep pumping until the pressure is stable. If the pressure does *not* stabilize, look for obvious leaks before you submerge the transducer in the water. The purpose of submerging the transducer is to find small leaks that are not otherwise detectable.
- 4** Once 150 mm Hg is reached, observe the pressure for 45 seconds.

| | |
|---|---|
|  Do not immerse if pressure drops | <p>Caution Test-c3</p> <p>If the pressure drops to zero after you use the pump, do not place the transducer in the tank.</p> |
|---|---|

| | |
|--|---|
|  Keep plug dry | <p>WARNING T-w5</p> <p>To prevent electrical shock and damage to the transducer, the connector pins in the transducer plug must always be completely dry before you connect to a system.</p> |
|--|---|

- 5 If the pressure appears stable, place the transducer and covered plug in the water tank for 45 seconds. Water temperature should be between +10 °C and +40 °C.
- 6 With the transducer in the water tank, observe carefully to make sure that no bubbles escape from the transducer.

NOTE: *Some small, insignificant bubbles may escape from the housing/gasket interface when the transducer is placed in the tank - this is not a leak. A leak is indicated by a constant, steady stream of escaping bubbles.*

| | |
|--|---|
|  If you see bubbles, do not release pressure | <p>Caution Test-c5</p> <p>If you see any bubbles, remove the transducer from the tank before you release the pressure.</p> |
|--|---|


- 7 If you find a leak, contact your BK Medical service representative to have the transducer repaired.
- 8 When leakage testing is finished, dry the transducer and connector before releasing the pressure and removing the test lid.
- 9 If the transducer will be further processed in liquid (for example, cleaned manually), attach the watertight plug lid (not the test lid) properly.

Leakage Testing Kits

There are two leakage testing kits, UA1404 and UA1414.

Check your product's Product Data Sheet for the appropriate leakage testing kit.

It is not likely that the tester will need cleaning, but you can wipe the lid with a mild detergent, then wipe it with tap water, and then dry with a soft cloth. Do not try to clean the pump.

| | |
|---|---|
|  Do not use test lid for reprocessing | <p>Caution Test-c2</p> <p>The lid part of the leakage tester is for testing only. Do not use it when you reprocess the transducer.</p> |
|---|---|

Yearly Preventive Maintenance and Performance Test

To ensure proper performance of the entire ultrasound system, preventive maintenance of the system, including a performance test, should be carried out once a year by a BK Medical technician or a suitably qualified engineer.

Follow local procedures or consult your BK Medical service representative about how to perform this check.


Circles and shadows when imaging in air

When you observe the monitor image from an array transducer that is not contacting a surface, you may see circles (or lines) and shadows. The concentric circles (or lines) are caused by re-reflections within the transducer and may not be uniform; they disappear when you image tissue. The shadows are caused by variations in the


transducer elements and the structure of the transducer surface. They do not indicate that the transducer is beginning to fail, and they do not influence general image performance.

All the equipment necessary for carrying out system preventive maintenance can be obtained from BK Medical.

Yearly Check of Type BF Transducers

| | |
|---|---|
|  Check of Type BF transducers | WARNING Check-w3 To prevent electrical shock, all transducers that comply with Safety Standard EN60601-1 (IEC60601-1) Type BF must be checked once a year to ensure that they still comply with the requirements of this standard. Transducers that need to be checked have the letters BF or the symbol printed on them. This check must be carried out by qualified personnel. Contact your BK service representative if you need any help checking your transducers. |
|---|---|

Chapter 6: Watertight Protection Devices

| | |
|---|---|
| <p style="text-align: center;">  Examine plug and waterproof protection for damage </p> | <p>Caution Plug-c4</p> <p>Before you reprocess the transducer, inspect the watertight protection device and the transducer plug. If you find any signs of damage, do not immerse the plug. If liquid comes into contact with the plug connector pins, the transducer may be destroyed.</p> <p>Examine the edges of the plug case that contact the lid and also the watertight protection device for cracks and marks. Examine the silicone seal of the plug lid. Look for deep scratches and grooves, holes or tears, brittleness, and looseness anywhere.</p> <p>The transducer or watertight protection device must be checked by a BK service representative if you find signs of damage.</p> |
|---|---|

Watertight Plug Lids

Some BK Medical transducers are supplied with a watertight plug lid that is designed to protect the transducer connector plug during immersion in liquids (including Medivators and STERIS system processing). When the lid is fitted, the entire transducer (including the covered plug) can be reprocessed using approved liquid-based methods.

There are three types of watertight plug lids: Type 1, Type 2 and the Watertight Plug Cover (see page 58).




Lid off



Lid on

Figure 3. Watertight Plug Lid Type 1.



Lid off



Lid on

Figure 4. Watertight Plug Lid Type 2.

2 covers
required for
8838

How to Attach and Detach Watertight Plug Lid Type 1

NOTE: The 8838 transducer has two plugs, one with a watertight plug lid, and a small one with a watertight plug cover. Both plugs must be protected correctly before the transducer is immersed.

To attach the watertight plug lid:

- 1 Place the lid on top of the plug with the locking screws outside. Make sure that the lid is properly aligned and seated on top of the plug.
- 2 Simultaneously turn both locking screws clockwise to attach the lid tightly to the plug. (You can also alternately tighten one screw and then the other, but do not fully tighten one before tightening the other). Tighten both screws until they cannot be tightened any more, but do not over-tighten them.

Do not fully tighten one screw before tightening the other screw, and do not put the lid on crooked, as these actions can damage the threads in the screw holes.

To detach the watertight plug lid:

- 1 Unscrew both locking screws by turning them counterclockwise.
- 2 Remove the lid from the plug.

How to Attach and Detach Watertight Plug Lid Type 2

To attach the watertight plug lid:

- 1 Place the lid on top of the plug with the locking pins unlocked (pointing at the open lock sign). Make sure that the lid is properly aligned and seated on top of the plug.
- 2 Turn both locking pins 90° clockwise to the locked position.

To detach the watertight plug lid:

- 1 Turn the locking pins 90° counterclockwise so they are unlocked.
- 2 Remove the lid from the plug.





Figure 5. Locking pins in locked and unlocked positions.

How to Attach the Watertight Plug Cover

Single-element transducers Type 2052 and 8838 come with a special watertight plug cover.

Protect plug
before
immersing

Screw the cover on tightly before you immerse the plug.



Lid off



Lid on

Figure 6. Watertight Plug Cover.

Chapter 7: Disposal


When you dispose of ultrasound equipment, you must follow national rules for the various materials in the equipment. Within the EU, you must send it to appropriate facilities for recovery and recycling.

BK Medical systems and transducers contain many different materials, but none require any special treatment compared with what would normally be expected for materials used in electronic equipment.

Be aware, however, that the printed circuit boards in the system are made of epoxy, the monitor's flat panel contains heavy metals and the system contains a small lithium battery.

For further information about the material composition of BK Medical equipment, contact your BK Medical service representative.

In general, dispose of the equipment in a way that minimizes the effects on the environment.

| | |
|--|--|
|  Contami- nated items | <p>WARNING D-w1</p> <p>For disposal of contaminated items such as transducer covers or needle guides or other disposable items, follow disposal control policies established for your office, department or hospital.</p> |
|--|--|

Packaging Material The packaging does not contain heavy metals or other dangerous materials. Follow your local procedure for disposing of and recycling non-dangerous waste.


Appendix: Reprocessing Information and Tables

Approved Disinfectants

In Germany and the USA, you must use cleaning and disinfection products that have been approved by the DGHM (Germany) or EPA or FDA (USA).

Transducer Compatibility

The **Transducer Reprocessing Methods** table shows which reprocessing methods your transducer can withstand, assuming that you follow the process manufacturer's instructions.

| | |
|---|--|
|  | <p>Caution: T-c5</p> <p>Using a non-recommended disinfection solution, an incorrect solution strength or immersing a transducer longer than recommended by the disinfectant manufacturer can damage the transducer.</p> |
|---|--|

NOTE: *Unless otherwise indicated, the table shows chemical (and physical) compatibility only – see the process manufacturer's own user instructions for information about the biological effectiveness of the method.*

Legend to the Transducer Reprocessing Methods table

| | |
|---------|---|
| • | means the transducer can withstand the process (when used according to manufacturer's instructions). |
| (blank) | means the transducer cannot withstand the process (or that it has not yet been tested with the process) |

Product and Process Manufacturers

| Product | Manufacturer |
|---|---------------------------------------|
| 3E-Zyme | Medisafe UK limited |
| Accel Prevention | Virox Technologies Inc. |
| Adaspor Single Shot | Cantel Medica (Italy) S.R.L. |
| Aniosyme DD1 | Laboratoires Anios |
| Antigermix S1 | Germitec |
| Astra VR | CIVCO Medical Solutions |
| Bodedex forte | BODE Chemie GmbH |
| Bomix Plus | BODE Chemie GmbH |
| Cavi Wipes/ CaviCide | Metrex Research |
| Cidex ADS/ Cidex OPA | Advanced Sterilization Products (ASP) |
| CIDEZYME/Enzol | Advanced Sterilization Products (ASP) |
| CIDEZYME XTRAMulti-Enzymatic Detergent/ CIDEZYME GL Enzymatic Detergent | Advanced Sterilization Products (ASP) |
| Cleanisept Wipes Forte | Dr. Schumacher GmbH |
| Clinell Sporidical wipes | GAMA Healthcare Ltd |
| Distel Wipes | Tristel Solutions Limited |
| Empower | Metrex |
| Ethanol 70% | None specific |
| Gigasept AF, Gigasept FF/Gigasept FF (Neu), Gigasept PAA concentrate | Schülke & Mayr GmbH |
| Gigazyme | Schülke & Mayr GmbH |
| Glutaraldehyde 2% - 3.4% | None specific |
| Incidin OxyFoam S Incidin OxyWipe S | Ecolab Deutschland GmbH |
| Intercept Wipes | MEDIVATORS Inc |
| Isopropanol 70% | MEDIVATORS Inc |
| Korsolex Basic | BODE Chemie GmbH |
| Korsolex Endo-Cleaner 0.5%/ Korsorex Endo-Disinfectant 1% | BODE Chemie GmbH |
| Korsolex Extra/ Korsorex Endo-Cleaner 0.5% | BODE Chemie GmbH |
| Metrizyme | Metrex |
| Mikrozoid AF liquid | Schülke & Mayr GmbH |
| Meliseptol Wipes sensitive | B. Braun Medical AG |
| Metricide/ Metricide 28/ Metricide Plus 30/ Metricide OPA Plus | Metrex |

| Product | Manufacturer |
|--|--|
| Mikrozoid PAA wipes | Schülke & Mayr GmbH |
| MATRIX Biofilm Remover | Whiteley |
| Meliseptol Foam | B. Braun Medical AG |
| Mikrobac Tissues | BODE Chemie GmbH |
| Neodisher Endo SEPT PAC/ Neodisher Septo PAC/ Neodisher Septo 3000/ Neodisher Endo Sept GA/ Neodisher MediClean Forte/ Neodisher Septo DN | Chemische Fabrik Dr. Weigert GmbH & Co. KG |
| Nu-Cidex | Advanced Sterilization Products (ASP) |
| OPAL | Whiteley Medical |
| Polystica 2x concentration. Enzymatic | STERIS Corporation |
| Perasafe/ Rely+On Perasafe | DuPont |
| Rapicide/ PA Disinfectant | MEDIVATORS Inc |
| Rapicide OPA/28 | MEDIVATORS Inc |
| RENO sterilizers (RENO-20, RENO-30, RENO-D50) - Cycle: Eco | Renosem Co., Ltd |
| RENO sterilizers (model RENO-S90, RENO-S130, RENO-S130D) -Cycle: Non-lumen for devices without hole/lumen - Cycle: Eco for devices with hole/lumen, e.g. 8666-RF, 9066 | Renosem Co., Ltd |
| RENO sterilizers (model RENO-S90, RENO-S130, RENO-S130D) -Cycle: Non-lumen for devices without hole/lumen - Cycle: Eco for devices with hole/lumen, e.g. 8666-RF, 9066 | Renosem Co., Ltd |
| Revital-Ox Resert/Resert XL HLD | STERIS Corporation |
| Revital-Ox detergent products: - Revital-Ox Bedside Complete -Revital-Ox 2X Concentrate Enzymatic Detergent - Revital-Ox Enzymatic Detergents | STERIS Corporation |
| Sani Cloth Plus Wipes/ Sani Cloth Super Wipes/ Sani-Cloth AF3 wipes | PDI |
| Sekusept Aktiv/ Sekusept MultiEnzyme P | Ecolab Deutschland GmbH |
| Steranios 2%, 2% N.G., 2% E.C.S | Laboratoires Anios |
| STERIS System 1/STERIS System 1E | STERIS Corporation |
| STERIS V-PRO 1/STERIS V-PRO 1 Plus/STERIS V-PRO 60/STERIS V-PRO maX | STERIS Corporation |
| STERIZONE VP4 | TSO ₃ Inc. |
| STERRAD 100S/ 200/ STERRAD NX/ STERRAD 100NX | Advanced Sterilization Products (ASP) |
| Suma Med Enzyme | Diversey |
| TD100 & TD5 | CS Medical |
| Thermosept PAA | Schülke & Mayr GmbH |
| Tristel Duo for Ultrasound/ Tristel Fuse for Instruments (Tristel Fuse for Stella)/ Tristel Trio Wipe System | Tristel Solutions Limited |
| Trophon EPR /Trophon2 | Nanosonics Ltd. |
| Wavicide 01 | Medical Chemical Corp. |
| Matachana 130LF, Webeco FA90, Webeco FA95 | Matachana Group |

Reprocessing Methods for flex Focus/bk3000/bk3500/bk5000/bkSpecto Systems

| BK Medical Systems - EXCEPT SCREEN . For bk3000, bk5000, bk3500, bkSpecto, and Flex Focus systems | | |
|--|-------------------------|-----------------------|
| Validated Products | Cleaning/Disinfection | Manufacturer |
| Tristel Duo for Ultrasound Wipes | Cleaning & Disinfection | Tristel Solutions Ltd |
| Ethanol 70% (wiping) | Disinfection | None specific |
| Isopropanol 70% (wiping) | Disinfection | None specific |
| Material Compatible Products ^a | Cleaning/Disinfection | Manufacturer |
| Any mild, non-corrosive, pH neutral (pH≈7) detergent product intended for wiping ultrasound devices/systems are considered to be material compatible with the system | Cleaning | None specific |
| Distel Wipes | Cleaning & Disinfection | Tristel Solutions Ltd |

a. Always follow product manufacturer's instructions, and use a slightly damp, soft, lint-free cloth or wipe to avoid liquid intrusion into the screen, keyboard or system. Do not pour liquid directly onto any part of the system.


| BK Medical Systems. For GLASS SCREENS on bk3000, bk3500, bk5000, Flex Focus 700, and Flex Focus 800 systems | | |
|--|-------------------------|-----------------------|
| Validated Products ^{a b} | Cleaning/Disinfection | Manufacturer |
| Tristel Duo for Ultrasound Wipes | Cleaning & Disinfection | Tristel Solutions Ltd |
| Ethanol 70% (wiping) | Disinfection | None specific |
| Isopropanol 70% (wiping) | Disinfection | None specific |
| Material Compatible Products ^a | Cleaning/Disinfection | Manufacturer |
| Any mild, non-corrosive, pH neutral (pH≈7) detergent product intended for wiping ultrasound devices/systems are considered to be material compatible with the system | Cleaning | None specific |
| Distel Wipes | Cleaning & Disinfection | Tristel Solutions Ltd |
| BK Medical Systems. For LCD SCREENS on bkSpecto, Flex Focus 200, Flex Focus 400, and Flex Focus 500 systems | | |
| Validated Products ^{a b} | Cleaning/Disinfection | Manufacturer |
| Ethanol 70% (minimum) (wiping) | Disinfection | None specific |

a. Always follow product manufacturer's instructions, and use a slightly damp, soft, lint-free cloth or wipe to avoid liquid intrusion into the screen. Do not pour liquid directly onto the screen.

b. Remove any remaining marks on the screen by wiping with Ethanol, and drying with a dry, lint-free, soft cloth while the screen is still moist.

Recommended Reprocessing Method for LCD Screens

- 1 Use a lint-free, soft cloth/wipe to wipe the screen (e.g. clean room wipe class: Level 100/ISO 5)
- 2 Wipe the screen in the direction shown in the figure - from the outside towards the center of the screen
- 3 If necessary use a soft, lint-free cloth/wipe, slightly moistened with an approved cleaner/disinfectant and wipe the screen
- 4 Remove any remaining marks on the screen by wiping with Ethanol, and drying with a dry, lint-free, soft cloth while the screen is still moist







LCD screen wiping direction

NOTE: Avoid liquid intruding under the screen.

| | |
|--------------------------|---|
| Screen contamination | <p>WARNING Reproc-w5</p> <p>To avoid contamination of the speaker area, do not touch this area when you turn the monitor to the horizontal or vertical position.</p> |
|--------------------------|---|











Validated Reprocessing Methods and Material Compatible Reprocessing Methods for the T7P2m (9027) Transducer

| Tip until 100 cm marker - submersible | Handle and cable part - not submersible, wiping only | Connector (and white flex relief) - not submersible, wiping only |
|---|---|---|
|  |  |  |
| Validated methods | | |
| 3E-Zyme | | |
| Revital-Ox Resert XL HDL | | |
| Material Compatible Methods | | |
| Mild pH neutral detergent | Mild pH neutral detergent wipe | Mild pH neutral detergent wipe |
| Aniosyme DD1 | Clinell Sporicidal wipes | Tristel Trio Wipes System |
| Cidex ADS, OPA and Nu-Cidex | Ethanol (EtOH) 70% wipes | Ethanol 70% wipes |
| Cidezyme / Enzol | Isopropyl alcohol (IPA) 70% wipes | IPA 70% wipes |
| Empower | Mikrozyd® PAA wipes | |
| Gigasept AF, Gigasept FF (Neu) and Gigasept PAA concentrate | Sani-Cloth AF3 wipes | |
| Korsorex extra | Tristel Trio Wipes System | |
| Metricide, Metricide 28, Metricide OPA Plus and Metricide Plus 30 | | |
| Metrizyme | | |
| Neodisher MediClean forte/ Neodisher Septo 3000 | | |
| Perasafe | | |
| Prolystica 2x conc. | | |
| Revital-Ox Resert XL HDL | | |
| Sekusept Aktiv | | |
| Steranios 2%, 2% N.G., 2% E.C.S | | |
| TD100 & TD5 | | |
| Tristel Trio Wipe System | | |
| Wavicide 01 | | |

| | |
|--|---|
|  Reprocessing | <p>WARNING Reproc-w4</p> <p>Extra care must be taken when cleaning this transducer, as there is no automatic cleaning system that can reprocess the entire transducer.</p> |
|--|---|

Material Compatible Reprocessing Methods for bkFusion Hardware

Do not immerse electronics unit, transmitter, sensor or cables.

| | | | | | |
|---|--|--|---|---|---|
|  | <p>EM Transmitter Stand (wheelbase and pole)</p> | <p>See "Reprocessing Methods for flex Focus/bk3000/bk3500/bk5000/bkSpecto Systems" on page 64</p> |  | <p>Mount for EM Transmitter</p> | <p>pH neutral (pH 6-8), non-corrosive cleaning products intended for wiping medical devices</p> <p>Ethanol wiping</p> |
|  | <p>Shelf for EM Control Unit</p> | <p>See "Reprocessing Methods for flex Focus/bk3000/bk3500/bk5000/bkSpecto Systems" on page 64</p> |  | <p>EM Transmitter Stand (articulated arm)</p> | <p>Cleaning:</p> <ol style="list-style-type: none"> 1. Wipe with a soft, non-abrasive cloth moistened with a mild, pH neutral (pH 6-8), non-abrasive detergent solution - or use a wipe manufactured for this purpose. 2. Wipe with soft cloth moistened with cold water 3. Wipe dry with a lint-free cloth <p>Disinfection:</p> <ol style="list-style-type: none"> 1. Use ethanol or isopropanol 70% wipes and spray |
|  | <p>EM Control Unit</p> | <p>See "Material Compatible Reprocessing Methods for Guides, Attachments and bkFusion Hardware" on page 73</p> | | | |
|  | <p>EM Sensor and Cable UA2371</p> | | | | |
|  | <p>EM Transmitter</p> | | | | |
|  | <p>Universal Bedside Clamp</p> | | | | |
|  | <p>Sensor Clamp UA2377</p> | | | | |
|  | <p>Sensor Clamp UA2399</p> | | | | |

Validated Reprocessing Methods for the bk3000/bk3500/bk5000/bkSpecto Transducer Series

| Validated Reprocessing Methods ^{ab} <i>(Follow local regulations for minimum reprocessing. Check table 4 on page 29)</i> | | Transducers | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--|--|--------------------------|------------|-------------|------------|-------------|---------------------------|--------------|--------------|-------------|--------------|------------|---------------|--------------|----------------|---------------|--------------|--------------|-------------|----------------|---------------|---------------|---------------|--------------|---------------|--------------|--------------------------|---------------|--------------|--------------|
| | | Surface | | | | | | | | | | | | Endocavity | | | | | | Intraoperative | | | | Neurology | | Phased Array | | Special | | |
| | | 5C1e ^c (9085) | 6C2 (9040) | 6C2s (9023) | 9C2 (9002) | 14L3 (9051) | 14L3e ^c (9086) | 13L4w (9011) | 10L2w (9022) | 18L5 (9070) | 18L5s (9081) | 8L2 (9032) | E14C4t (9018) | E14C4 (9067) | E14CL4b (9048) | E11C3b (9008) | E13C2 (9029) | E10C4 (9019) | 20R3 (9052) | I14C5I (9015) | I14C5T (9016) | I12C5b (9024) | I12C4f (9066) | N13C5 (9062) | N11C5s (9063) | 5P1 (9077) | 5P1e ^c (9087) | X18L5s (9009) | X14L4 (9038) | X12C4 (9026) |
| Manual Cleaning | 3E-Zyme | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . |
| | Manual Disinfection | Korsorex Basic | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . |
| Ethanol 70% (wiping) | | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| Revital-Ox Resert/Resert XL HLD | | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| Tristel Fuse For Instruments (Tristel Fuse for Stella) | | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| Automated Disinfection | Medivators Advantage Plus Intercept (detergent), Rapicide PA Disinfectant flush: 70% isopropyl alcohol | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | 2 | . | . | |
| Sterilization | STERIS System 1, 1 Plus and 1 Express ^d STERIS System 1E | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | STERIS V-Pro 1 Plus, V-Pro 60 Non lumen cycle | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | STERIS V-Pro maX Non lumen cycle or Flexible cycle | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Sterrad 100NX Standard cycle | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Sterrad 100NX Express cycle | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Sterrad NX Standard cycle | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Sterrad 100S One cycle only (USA), Short cycle (rest of the world) | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Sterrad 200 Short cycle | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| Matachana 130LF, Webeco FA90, Webeco FA95 (60°C cycle) | | | | | | | | | | | 1 | 1 | 1 | | 1 | 1 | | 1 | 1 | 1 | | | | | | | | | | |

a. Cleaning is an important first step for all transducer types. Every transducer used in a critical procedure must be sterilized as the final reprocessing step.

b. Follow product manufacturer's instructions and do not exceed transducer-specified limits. See Table 3 on page 25 for more information.

c. Transducer connector is not immersible.

d. STERIS SYSTEM 1, SYSTEM 1 Plus and SYSTEM 1 Express have not been market cleared by the FDA.

1. Use a suitable connector lid and valve from the system manufacturer. See page 44.

2. Use of Medivators Advantage Plus with X18L5s has not been CE-marked.

| Material Compatible Reprocessing Methods ^{abc} (Follow local regulations for minimum reprocessing. Check table 4 on page 29) | | Transducers | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--|--|--------------------------|------------|-------------|------------|-------------|---------------------------|--------------|--------------|-------------|--------------|------------|---------------|--------------|----------------|---------------|--------------|----------------|-------------|---------------|---------------|---------------|---------------|--------------|---------------|------------|--------------------------|---------------|--------------|--------------|---|---|---|---|---|---|---|
| | | Surface | | | | | | | | | | Endocavity | | | | | | Intraoperative | | | Neurology | | Phased Array | | Special | | | | | | | | | | | | |
| | | 5C1e ^d (9085) | 6C2 (9040) | 6C2s (9023) | 9C2 (9002) | 14L3 (9051) | 14L3a ^e (9086) | 13L4w (9011) | 10L2w (9022) | 18L5 (9070) | 18L5s (9081) | 8L2 (9032) | E14C4t (9018) | E14C4 (9067) | E14CL4b (9048) | E11C3b (9008) | E13C2 (9029) | E10C4 (9019) | 20B3 (9052) | 114C5i (9015) | 114C5T (9016) | 112C5b (9024) | 112C4f (9066) | N13C5 (9062) | N11C5s (9063) | 5P1 (9077) | 5P1e ^d (9087) | X18L5s (9009) | X14L4 (9038) | X12C4 (9026) | | | | | | | |
| Sterilization | OPAL | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Rapicide/Rapicide OPA/28 | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Rely+On Perasafe | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Revital-Ox Resert / Resert XL HLD (Automated use) | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Sani Cloth Plus Wipes/Sani Cloth Super Wipes | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Sekusept Aktiv | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Steranios 2%, 2% N.G., 2% E.C.S | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Thermosept PAA | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . |
| | Tristel Duo for Ultrasound/Tristel Trio Wipe System | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . |
| | Trophon EPR / Trophon2 (Fit transducer completely above line mark) | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . |
| Sterilization | RENO-20/RENO-30/RENO-D50: ECO Cycle | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | RENO-S90, RENO-S130, RENO-S130D: Non Lumen Cycle | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . |
| | RENO-S90, RENO-S130, RENO-S130D: ECO Cycle | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . |

- a. Cleaning is an important first step for all transducer types. Every transducer used in a critical procedure must be sterilized as the final reprocessing step.
- b. Follow product manufacturer’s instructions and do not exceed transducer-specified limits. See Table 3 on page 25 for more information.
- c. This table of reprocessing methods only indicates that BK Medical has evaluated the transducers for material compatibility, when reprocessed according to the product/system IFU, and not for their efficacy in attaining the appropriate level of disinfection. Ensure that the reprocessing chamber is an appropriate size for the transducer.
- d. Transducer connector is not immersible.

Validated Reprocessing Methods for the Flex Focus and Pro Focus Transducer Series

| Validated Reprocessing Methods ^{ab} <i>(Follow local regulations for minimum reprocessing. Check table 4 on page 29)</i> | | Transducers | | | | | | | | | | | | | | | | | | | | | | | | |
|--|---|----------------|------|------|------|------|------|------|------|------|------------|---------|------|-------|------|-------------------|------|------|-------------------|------|-------|------|------|-------------------|-------------------|------|
| | | Intraoperative | | | | | | | | | Endocavity | | | | | Surface | | | | | | | | | | |
| | | 8666-RF | 8809 | 8815 | 8816 | 8824 | 8826 | 8836 | 8862 | 8863 | 2052 | 8667 | 8808 | 8808e | 8818 | 8819 ^c | 8838 | 8848 | 8670 ^c | 8811 | 8820e | 8822 | 8823 | 8830 ^c | 8837 ^c | 8870 |
| Manual Cleaning | 3E-Zyme | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . |
| | Manual Disinfection | Korsolex Basic | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . |
| Ethanol 70% (wiping) | | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| Revital-Ox Resert/Resert XL HLD | | | | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| Tristel Fuse for Instruments (Tristel Fuse for Stella) | | | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| Automated Disinfection | Medivators Advantage Plus Intercept (detergent), Rapicide PA Disinfectant, flush: 70% isopropyl alcohol | | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| Sterilization | STERIS System 1,1 Plus, and 1 Express ^d STERIS System 1E | . | . | . | . | . | . | . | 4 | 4 | . | . | 5 | . | . | . | . | . | . | . | . | . | . | . | . | |
| | STERIS V-Pro 1 Plus, V-Pro 60 Non lumen cycle | . | 2 | . | . | . | . | . | . | . | . | 9 | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | STERIS V-Pro maX Non lumen cycle or Flexible cycle | . | 2 | . | . | . | . | . | . | . | . | 9 | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | STERIZONE VP4 | 1 | . | . | . | . | . | . | . | . | . | . | . | 6 | . | . | 7 | . | . | . | . | . | . | . | . | |
| | Sterrad NX & 100NX Standard cycle | 1 | 2 | 3 | 3 | 3 | . | . | . | . | . | . | . | 6 | . | . | 7 | . | . | . | . | . | . | . | . | |
| | Sterrad 100S One cycle only (USA), Short cycle (rest of the world) | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Sterrad 200 Short cycle | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | |
| | Matachana 130LF, Webeco FA90, Webeco FA95 (60°C cycle) | | | 8 | 8 | 8 | . | . | . | . | . | 8/ 9 | 8 | 8 | . | . | 8 | . | . | . | . | . | . | . | . | |

a. Cleaning is an important first step for all transducer types. Every transducer used in a critical procedure must be sterilized as the final reprocessing step.

b. Follow product manufacturer's instructions and do not exceed transducer-specified limits. See Table 3 on page 25 for more information.

c. Transducer connector is not immersible.

d. STERIS SYSTEM 1, SYSTEM 1 Plus and SYSTEM 1 Express are not market cleared by the FDA.

1. Valid for transducers with a serial number higher than 1911238. If in doubt, contact your BK representative.

2. Valid for transducers with a serial number higher than 3991000. If in doubt, contact your BK representative.

3. Valid for transducers with a serial number higher than 1910000. If in doubt, contact your BK representative.

4. In Canada, do not use STERIS SYSTEM 1E for this transducer.

5. Use together with Quick Connect QMC1733INT and tray: C1220S1 for STERIS SYSTEM 1; Quick Connect QMC1733INT and tray: C1220INT for SYSTEM 1 Plus / 1 Express and Quick Connect QMC1733E and Tray: C1220E with SYSTEM 1E.

6. Valid for transducers with a serial number higher than 3890001. If in doubt, contact your BK representative.

7. Valid for transducers with a serial number higher than 3900001. If in doubt, contact your BK representative.

8. Use a suitable connector lid and valve from the system manufacturer. See page 44.

9. Valid for transducers with a serial number higher than 3018205. If in doubt, contact your BK representative.

Material Compatible Reprocessing Methods for the Flex Focus and Pro Focus Transducer Series

| Material Compatible Reprocessing Methods ^{abc} (Follow local regulations for minimum reprocessing. Check table 4 on page 29) | | Transducers | | | | | | | | | | | | | | | | | | | | | | | | |
|--|---|----------------|------|------|------|------|------|------|------|------------|------|------|------|-------|------|-------------------|---------|------|-------------------|------|-------|------|------|-------------------|-------------------|------|
| | | Intraoperative | | | | | | | | Endocavity | | | | | | | Surface | | | | | | | | | |
| | | 8666-RF | 8809 | 8815 | 8816 | 8824 | 8826 | 8836 | 8862 | 8863 | 2052 | 8667 | 8808 | 8808e | 8818 | 8819 ^d | 8838 | 8848 | 8670 ^d | 8811 | 8820e | 8822 | 8823 | 8830 ^d | 8837 ^d | 8870 |
| Cleaning | 3E-Zyme (Automated use) | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Bodedex forte | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | CIDEZYME XTRA Multi-Enzymatic Detergent/ CIDEZYME GL Enzymatic Detergent | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | pH neutral (pH 6-8), non-corrosive cleaning products intended for medical devices | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Intercept Wipes/Intercept Detergent | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Gigazyme | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Korsolex Endo-Cleaner 0,5% | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | MATRIX Biofilm Remover | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Neodisher MediClean Forte | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Prolystica 2x concentration. Enzymatic | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Revital-Ox Bedside Complete/2X Concentrate Enzymatic Detergent/Enzymatic Detergents | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Sekusept MultiEnzyme P | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| Suma Med Enzyme | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | |
| Disinfection (continued on next page) | Accel Prevention (wipes, ready-to-use liquid, concentrate) | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Adaspor Single Shot | | | | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Antigemix S1 | | | | | | | | | | | | * | * | * | * | * | * | * | * | * | * | * | * | * | |
| | Astra VR (with approved disinfectant) | | | | | | | | | | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Bomix Plus | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Cavi Wipes/CaviCide | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Cidex OPA | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Cleanisept Wipes Forte | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Isopropanol 70% (wiping) | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Gigasept FF | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Glutaraldehyde 2% – 3.4% | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Incidin OxyFoam/Incidin OxyWipe S | | | | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Korsolex Endo Disinfectant 1%/Korsolex Extra | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Korsolex Basic (Automated use) | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Meliseptol Foam | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Meliseptol Wipes Sensitive | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Metricide OPA Plus | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Mikrobac Tissues | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Mikrozoid AF Liquid | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| | Neodisher Endo Sept GA | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * |
| Neodisher Septo DN | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | |
| Nu-Cidex | | | | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | * | |
| Table continues on next page | | | | | | | | | | | | | | | | | | | | | | | | | | |

- a. Cleaning is an important first step for all transducer types. Every transducer used in a critical procedure must be sterilized as the final reprocessing step.
- b. Follow product manufacturer’s instructions and do not exceed transducer-specified limits. See Table 3 on page 25 for more information.
- c. This table of reprocessing methods only indicates that BK Medical has evaluated the transducers for material compatibility, when reprocessed according to the product/system IFU, and not for their efficacy in attaining the appropriate level of reprocessing. Ensure that the reprocessing chamber is an appropriate size for the transducer.
- d. Transducer connector is not immersible.

1. Transducers with a serial number higher than 1912156 can be reprocessed with STERIS V-Pro models. In in doubt, contact your BK representative.

| Material Compatible Reprocessing Methods ^{abc} <i>(Follow local regulations for minimum reprocessing. Check table 4 on page 29)</i> | | Transducers | | | | | | | | | | | | | | | | | | | | | | | | |
|---|---|----------------|------|------|------|------|------|------|------|------|------------|------|------|-------|------|-------------------|---------|------|-------------------|------|-------|------|------|-------------------|-------------------|------|
| | | Intraoperative | | | | | | | | | Endocavity | | | | | | Surface | | | | | | | | | |
| | | 8666-RF | 8809 | 8815 | 8816 | 8824 | 8826 | 8836 | 8862 | 8863 | 2052 | 8667 | 8808 | 8808e | 8818 | 8819 ^d | 8838 | 8848 | 8670 ^d | 8811 | 8820e | 8822 | 8823 | 8830 ^d | 8837 ^d | 8870 |
| | OPAL | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Rapicide | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | RAPICIDE OPA/28 | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Rely+On Perasafe | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Revital-Ox Resert / Resert XL HLD (Automated use) | | | • | • | | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Sani Cloth Plus Wipes/Super Wipes | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Sekusept Aktiv | | | • | • | | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Steranios 2%, 2% N.G., 2% E.C.S | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Thermosept PAA | | | • | • | | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Tristel Duo for Ultrasound | | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Tristel Trio Wipe System | | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Trophon EPR / Trophon2 <i>(Fit transducer completely above line mark)</i> | | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| Sterilization | RENO-20/RENO-30/RENO-D50: ECO Cycle | • | | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | RENO-S90, RENO-S130, RENO-S130D: Non Lumen Cycle | | | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | RENO-S90, RENO-S130, RENO-S130D: ECO Cycle | • | | | | | | | | | | | | | | | | | | | | | | | | |

- a. Cleaning is an important first step for all transducer types. Every transducer used in a critical procedure must be sterilized as the final reprocessing step.
- b. Follow product manufacturer's instructions and do not exceed transducer-specified limits. See Table 3 on page 25 for more information.
- c. This table of reprocessing methods only indicates that BK Medical has evaluated the transducers for material compatibility, when reprocessed according to the product/system IFU, and not for their efficacy in attaining the appropriate level of reprocessing. Ensure that the reprocessing chamber is an appropriate size for the transducer.
- d. Transducer connector is not immersible.

1. Transducers with a serial number higher than 1912156 can be reprocessed with STERIS V-Pro models. In in doubt, contact your BK representative.

Material Compatible Reprocessing Methods for Guides, Attachments and bkFusion Hardware

| Material Compatible Reprocessing Methods ^a | Attachments, Needle Guides | | | | | | | | | | | | | | bkFusion Hardware | | | | | | |
|---|---|--------|--------|--------|--------|--------|--------|--------|--------|--------|-------------------|--------|--------|--------|---------------------|----------------------------------|-------------------------------------|-------------------------|-------------------------|--------------------------|---|
| | UA0671 | UA0672 | UA0673 | UA1232 | UA1239 | UA1250 | UA1251 | UA1256 | UA1282 | UA1324 | UA1325 / UA1325-W | UA1326 | UA1327 | UA1328 | UA2377 Sensor Clamp | UA2399 Sensor Clamp ^b | UA2371 EM Sensor and Cable (wiping) | Universal Bedside Clamp | EM Transmitter (wiping) | EM Control Unit (wiping) | |
| Cleaning ^c | pH neutral (pH 6-8), non-corrosive cleaning products intended for medical devices | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Intercept Detergent | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | • | | | | |
| | 3E-Zyme | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Gigazyme | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Korsolex Endo-Cleaner 0,5% | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Neodisher MediClean Forte | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Sekusept MultiEnzyme P | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| Disinfection ^c | Accel Prevention (wipes, ready-to-use liquid, concentrate) | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Adaspor Single Shot | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Bomix Plus | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Cidex OPA | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Gigasept FF | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Glutaraldehyde 2% - 3.4% | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Isopropanol 70% / Ethanol 70% (wiping) | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Korsolex Endo Disinfectant 1% | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Korsolex Extra | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Korsolex Basic | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Rapicide PA | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Metricide OPA Plus | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Neodisher Endo Sept GA | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Neodisher Endo SEPT PAC/Neodisher Septo PAC | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Neodisher Septo DN | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Nu-Cidex | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | OPAL | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Rapicide | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Rapicide OPA/28 | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Rely+On Perasafe | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| Revital-Ox Resert/Resert XL HLD | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | | |
| Sekusept Aktiv | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | | |
| Thermosept PAA | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | | |
| Tristel Fuse For Instruments | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | | |
| Sterilization ^c | Ethylene Oxide (ETO) | | | | | | | | | | | | | | | | | | | | • |
| | STERIS System 1, 1E, 1 Plus and 1 Express | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | STERIS V-Pro 1 Plus/ V-Pro 60/V-Pro maX | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| | Sterrad NX/100NX/100S/200 | • | • | • | | • | • | • | • | • | • | • | • | • | • | • | | | | | |
| Autoclaving (See "Autoclaving" on page 46) | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |

- a. This table of reprocessing methods only indicates that BK Medical has evaluated the guides and attachment for material compatibility and not for their efficacy in attaining the appropriate level of reprocessing.
- b. Separate the clamp into two components before reprocessing.
- c. Follow product manufacturer's instructions.

1. Transducers with a serial number higher than 1912156 can be reprocessed with STERIS V-Pro models. In in doubt, contact your BK representative.

Reprocessing Methods for Remote Controls UA1237 and UA2361

NOTE: Remote Control UA1237: Only remote controls with a serial number (S/N#) of 1000200 or higher can be immersed. For further information, see “Reprocessing Remote Controls” on page 36.

| | | Reprocessing Methods | Lid and Batteries |
|--|---------------------|---|---|
| Material Compatible Methods ^a | Manual Cleaning | Non-corrosive, pH neutral (pH 6-8) cleaning products intended for medical devices | Lid on |
| | Manual Disinfection | Ethanol 70% (wiping) | Batteries inside the remote control Lid on |
| | | Cavi Wipes | |
| | | CaviCide | |
| | | Cidex OPA | |
| | | Gigasept FF | |
| | | Glutaraldehyde 2% | |
| | | Isopropanol 70% (wiping) | |
| | | Korsolex Basic | |
| | | Korsolex Endo Disinfectant 1% | |
| | | Korsolex Extra | |
| | | Meliseptol Foam | |
| | | Metricide OPA Plus | |
| | | Nu-Cidex | |
| | | Rapicide OPA/28 | |
| | | Rely+On Perasafe | |
| Sani Cloth Plus Wipes | | | |
| Sani Cloth Super Wipes | | | |
| Validated methods ^{a,b,c} | Sterilization | STERIS SYSTEM 1 STERIS SYSTEM 1 Plus STERIS SYSTEM 1 Express STERIS SYSTEM 1E | Batteries inside the remote control Lid on |
| | | STERIS V-PRO 1 <i>Standard Cycle</i> | Batteries <i>with</i> but not <i>inside</i> the remote control Lid off |
| | | STERIS V-PRO 1 Plus STERIS V-PRO maX STERIS V-PRO 60 <i>Lumen cycle</i> | Batteries <i>with</i> but not <i>inside</i> the remote control Lid off |
| | | STERRAD 100S AND 200 <i>Standard Cycle (USA)</i> <i>Short Cycle (rest of world)</i> | Batteries <i>with</i> but not <i>inside</i> the remote control Lid off |
| | | STERRAD NX STERRAD 100NX <i>Standard Cycle</i> | Batteries <i>with</i> but not <i>inside</i> the remote control Lid off |

a. Follow product manufacturer's instructions and do not exceed remote control-specified limits.

b. STERIS SYSTEM 1, SYSTEM 1 PLUS, and SYSTEM1 Express are not market cleared by the FDA.

c. Because of material degeneration, STERRAD systems processing has been validated for a maximum of 50 reprocessing cycles with the remote control.

Overview of InstruSafe Instrument Protection Trays for BK Devices

The listed InstruSafe protection trays can be used for the storage, transportation and efficient sterilization of the majority of BK transducers with V-PRO and STERRAD sterilizers. Some V-PRO and STERRAD sterilizers have a small chamber, so ensure tray dimensions fit the processing chamber of your sterilizer.

Always follow manufacturers' instructions.

Prior to use, check transducer compatibility “Validated Reprocessing Methods for the bk3000/bk3500/bk5000/bkSpecto Transducer Series” on page 67, and “Material Compatible Reprocessing Methods for the bk3000/bk3500/bk5000/bkSpecto Transducer Series” on page 68, plus tray compatibility^a, with exact sterilization models and cycles.

| TRANSDUCER | TRAY ^a | DIMENSIONS |
|------------------------------------|-------------------|---|
| 20R3, 2052 | UA2431 | 63.18 cm (L) x 34.29 cm (W) x 6.03 cm (H) |
| I12C4f, 8666-RF | UA2432 | 63.18 cm (L) x 34.29 cm (W) x 6.03 cm (H) |
| X18L5s, 8809 | UA2433 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| I14C5T, 8816 | UA2434 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| I12C5b, 8824 | UA2435 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| X12C4, 8826 | UA2436 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| 9C2, 18L5, 13L4w, 8670, 8811, 8870 | UA2437 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| E14CL4b, E10C4, 8819, 8848 | UA2438 | 63.18 cm (L) x 34.29 cm (W) x 6.03 cm (H) |
| N13C5, 8862 | UA2439 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| 8836 | UA2440 | 60.64 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| N11C5s, 8863 | UA2441 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| 6C2s, 8823 | UA2442 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| E14C4t, 8808e, 8808, 8818 | UA2443 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| I14C5I, 8815 | UA2444 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| 5P1, 8837 | UA2445 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| 8L2, 14L3 | UA2446 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| E14C4, 8667 | UA2447 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| 6C2, 10L2w, 8820e, 8822, 8830 | UA2448 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |
| X14L4, 8838 | UA2449 | 49.85 cm (L) x 24.45 cm (W) x 6.03 cm (H) |

a. Manufactured by Summit Medical. Refer to www.instrusafe.com for further information, approved systems, cycles and trays.



Reprocessing Table for Craniotomy Transducer 8862 and Burr-Hole Transducer 8863

| | |
|------------------------|---|
| Product Name: | Craniotomy Transducer 8862 and Burr-Hole Transducer 8863 |
| Manufacturer: | BK Medical, Mileparken 34, 2730 Herlev, Denmark |
| Product Number: | Type 8862 and Type 8863 |
| Contact: | Your local BK Representative or info@bkmedical.com Tel. +45 4452 8100 |

This table tells you where to find specific information about reprocessing these transducers.

| | |
|---------------------|--|
| Abbreviations used: | C&C: <i>Care and Cleaning</i> |
| | UG: Transducer User Guide |

Reprocessing Information for Neurosurgical Transducers 8862 & 8863

| Process | Process Stage | Process Step | Recom'd | Not Recom'd | Specific information to be provided by manufacturer (attach details) | Document | Section |
|------------------------------------|--------------------------------|------------------------------------|---------|-------------|---|----------|--|
| Preparation at point of use | Washing after use | | X | | Specify type of detergent or agent to use for soak (for example [e.g.] alkaline, acidic, neutral pH, enzymatic solution, enzymatic foam, or water). | C&C | Pre Cleaning (Point-of use Processing). |
| | | Rinsing | | | Note: Soaking is not recommended. Rinse under running water. | C&C | Pre Cleaning (Point-of use Processing). |
| Decontamination | Preparation | Disassembly | | | Device specific disassembly instructions with pictures. | UG | Detailed diagrams in user guide show how needle guides click on and off. |
| | Cleaning (Includes rinsing) | Manual cleaning | X | | Specify any special cleaning brushes or tools needed. Specify water quality needed. Specify type of agent to use for cleaning (e.g. alkaline, acidic, neutral pH, enzymatic solution, enzymatic foam, or water). Specify minimum volume of water needed for rinsing. | C&C | Manual Cleaning. |
| | | Automated (Machine) Cleaning | | X | | | Automatic Reprocessing Methods. |
| | | Ultrasonic Cleaning | | X | | | |

(Sheet 1 of 3)

Reprocessing Information for Neurosurgical Transducers 8862 & 8863

| Process | Process Stage | Process Step | Recom'd | Not Recom'd | Specific information to be provided by manufacturer (attach details) | Document | Section |
|--------------------------------|-----------------|----------------|---------|-------------|---|----------|--|
| Disinfection | Liquid Chemical | Manual | X | | Specify compatible liquid chemicals that can be used. | C&C | Material Compatible Reprocessing Methods. |
| | | | X | | Specify validated exposure time to liquid chemical. | PD | Manual Disinfection. |
| | | | | | Specify water quality for rinse and minimum volume for rinsing. | C&C | Manual Disinfection. |
| Drying | Thermal | Automated Only | | X | | | |
| | | | X | | Specify how device should be dried (e.g. pressurized air at recommended maximum air pressure, manual wiping, heat, etc.). | C&C | Manual Disinfection. |
| Preparation and Packing | | | | | Specify maximum temperature the medical device can withstand. | | Automatic Reprocessing Methods. |
| | Reassembly | | | X | | | No reassembly before sterilization. |
| Maintenance | | | X | | Specify any requirements for ensuring functionality, e.g. sharpening, lubrication, testing device function, testing sheath integrity. | C&C | Checking and Maintaining Ultrasound Equipment. |
| Steam Sterilization | | | | X | | | |

(continued) (Sheet 2 of 3)

Reprocessing Information for Neurosurgical Transducers 8862 & 8863

| Process | Process Stage | Process Step | Recom'd | Not Recom'd | Specific information to be provided by manufacturer (attach details) | Document | Section |
|--|----------------------------------|--------------|---------|-------------|---|----------------|--|
| EO Sterilization | | | | X | | | |
| Other Sterilization Processes | | | | | Specify sterilization process including cycle and conditions for which device has been validated. | | Automatic Reprocessing Methods. |
| | STERRAD® | | X | | | C&C | STERRARD Systems. |
| | STERIS™ SYSTEM 1 Amsco™ V-PRO | | X X | | There is only one type of sterilization cycle. | C&C | STERIS SYSTEM 1 and 1E. STERIS Amsco V-PRO Systems. |
| Device to be Sterilized in Container Provided by Manufacturer | | | | X | | | |

(continued) (Sheet 3 of 3)

Index

Numerics

10L2w 67, 68, 69, 71
13L4w 67, 68, 69, 71, 73
14L3 67, 68, 69, 71, 73
18L5 67, 68, 69, 71
2050 58
2052 58
20R3 67, 68, 69
5P1 67, 68, 69, 71
6C2 67, 68, 69, 70, 71, 72, 73
6C2s 67, 68, 69, 70, 71, 72
8666-RF 49
8809 49
8838 57
9C2 67, 68, 69

A

Atmospheric limitation (symbol) 14
atmospheric pressure limits 24
autoclaving
 about 46
 excessive heat and non-steel puncture guides 49
 puncture attachments 46
 required for puncture attachments 34
 transducer accessories 46

B

batch code (symbol) 13
batch number (symbol) 13
biopsy channels, cleaning 33

C

Canada, licensed sheaths 25
caution (symbol) 12
checking
 for leaks 52
 overview 48
 performance test 53
 plug and cover before immersion 49
 transducer before immersion 50
 transducers and PCU 49
circles and shadows, imaging in air 53
Cleaning and Disinfecting the System 35
condensation 7, 24
consult instructions for use (symbol) 12
Covers and Gels 25
Creutzfeldt-Jakob disease 7, 28
critical devices, definition 29
cross-contamination 6

D

date of manufacture (symbol) 14
disinfection
 magnetic wheel mover 34
disinfection levels, table 29

disposal 12, 60
do not resterilize (symbol) 13
do not reuse (symbol) 13
dummy channel bracket 34
dust, avoid 7, 24

E

E10C4 67, 68, 69
E14C4 67, 68, 69, 70, 71, 73
E14C4t 67, 68, 69
E14CL4b 67, 68, 69
environmental limits
 system 24
 transducers 25
EO
 symbol for sterile 13
expiration date (symbol) 14

G

gas processing
 pressure 25

H

handle with care (symbol) 13
heat, avoid 7, 24
Humidity limitation (symbol) 14
humidity, storage 25

I

I12C4f 49
InstruSafe Instrument Protection Trays 75

L

leakage tester
 picture of test setup 52
leakage tester UA1404
 description 53
Leakage Testing Frequency Table 51
Leakage Testing Procedure 52
leaks, testing for 52
lot number (symbol) 13
lumens, cleaning 33

M

magnetic wheel mover
 cleaning 34
 disinfection 34
Manual Disinfection (Immersion or Wiping) 33
manufacturer (symbol) 12
markings, table of 15
Medivators® AER system 38

N

N13C5 67, 68, 69
noncritical devices, definition 29

not watertight (symbol) 15

O

obturator 34
operating environment
 system 24
 transducer 25

O-rings 34

P

Package damage (symbol) 14
packaging, disposing of 60
performance test 53
pressure limits
 gas processing 25
 operating 25
probe sheaths. See transducer covers.
puncture attachments
 autoclaving 46
 cleaning before sterilizing 30
 with sterile transducer cover 26
puncture guides
 non-steel and excessive heat 49

R

rectosigmoidoscope 34
REF (symbol) 13
remote control
 immersion (symbols) 15
Reprocessing Equipment 30
Reprocessing Other Accessories 34
Reprocessing the Remote Control 36

S

sealing (symbol) 13
semi-critical devices, definition 29
single-use (symbol) 13
steam sterilization. See autoclaving.
sterile (symbol) 13
STERILE EO (symbol) 13
STERIS Amsco V-PRO Systems 42, 43
STERIS systems 39, 40, 78
STERIZONE® VP4 43
STERRAD (symbol) 15
storage
 system, environmental limits 24
 transducer, environmental limits 25
 transducers 23
sunlight, avoid 7, 24
sunlight, avoid (symbol) 14
symbols, table of 15
system
 cleaning 35

T

temperature
 operating 25
 storage 25
Temperature limitation (symbol) 14

testing. See checking.

training

 before cleaning and disinfecting 28

transducer covers

 reducing cross-contamination risk 25

transducer holder

 cross-contamination 6

transducers

 storing 23

Type B (symbol) 13

Type BF (symbol) 13

Type BF transducers, yearly check 54

Type BF, defibrillator-proof (symbol) 13

U

UA 0513 34
UA0671 34
UA0672 34
UA0673 34
UA0674 34
UA1237 Remote Control 20
UA1237 Reprocessing Methods 74
UA1325-w 34
UA1404 50
UA1414 50, 52
UA1414 Leakage Testing Kit 50
UA2361 Remote Control 19
UA2361 Reprocessing Methods 74
UA2431 75
UA2432 75
UA2433 75
UA2434 75
UA2435 75
UA2436 75
UA2437 75
UA2438 75
UA2439 75
UA2440 75
UA2441 75
UA2442 75
UA2443 75
UA2444 75
UA2445 75
UA2446 75
UA2447 75
UA2448 75
UA2449 75
USA
 caution 7, 18
 market-cleared transducer covers 25
use by date (symbol) 14

W

warning (symbol) 12
waste disposal symbol 12
water standoff collar 34
Watertight Plug Cover 58
watertight plug cover 58
watertight plug lid for 90XX series transducers and on

| | |
|--|------------------------|
| 56 | WEEE waste (symbol) 12 |
| Watertight Plug Lids 56 | |
| Webeco FA90, Webeco FA95 and Matachana 130LF | X |
| 44 | X18L5s 49 |

 BK Medical ApS, Mileparken 34, 2730 Herlev, Denmark. T +45 4452 8100 F +45 4452 8199



BK Medical
8 Centennial Drive
Peabody
MA01960
USA
T + 1 978-326-1300
bkmedical.com

USA
Sales & Service
BK Medical
8 Centennial Drive
Peabody
MA01960
USA
T + 1 978-326-1300
F + 1 978-326-1399
bkmedical.com

Europe and Rest of World
Sales, Service & Design Center
BK Medical
Mileparken 34
2730 Herlev
Denmark
T +45 4452 8100
F +45 4452 8199
bkmedical.com