

Indications for use

Fetal Pillow is intended to elevate the fetal head and facilitate delivery of the fetus in women requiring a Caesarean Section at full dilation or those requiring a Caesarean Section after a failed instrumental vaginal delivery. Fetal Pillow is indicated for use in gestational age >37 weeks.

Caution

Federal (U.S.A) Law restricts the use of this device by or on order of a physician.

Warnings

The safety and effectiveness of Fetal Pillow has not been established in the following:

1. In women who have had a previous caesarean section
2. In women with a pregnancy less than 37weeks
3. Non vertex presentation
4. Pregnancy with Intra-uterine Fetal Death
5. Pregnancy induced hypertension
6. Intra-uterine Growth Retardation
7. Diabetes in pregnancy
8. Major congenital abnormalities
9. Presence of chorioamnionitis
10. Multiple gestations

Contraindications

Fetal Pillow should not be used in the presence of active genital infection, as it could increase the risk of ascending infection.

Precautions


1. DO NOT use air to inflate the device
2. Maximum inflation should not be more than 180cc
3. The device will fail to inflate if the dome/balloon surface of device is not in contact with fetal head when inserted
4. Make sure that the package is intact before use
5. Inflate the device with 60 cc saline prior to use to check the integrity of the device

Please read all information carefully

Failure to properly follow instructions may result in improper functioning of the device.

Device description

Fetal Pillow is a sterile single use device consisting of a base plate and a dome (inflatable balloon) made of silicone. A 100cm long tubing is attached to this for inflation. The tubing has a two-way tap at the distal end for inflation and deflation. A sterile 60cc syringe is provided with the device for inflation using sterile saline. The dome inflates only in upward direction when placed correctly.



Step 1:


Before Inserting

Insertion and inflation of the device should be carried out just before performing the Caesarean Section.

Inflate the device with 60 cc saline prior to use to check the integrity of the device. Empty the device using the syringe provided before insertion.

Hold the base plate of Fetal Pillow between fingers and thumb as shown and fold to squeeze the Dome (balloon) between the base-plate.

The tube attachment should be at the superior end during insertion as shown. If the tube attachment is facing downwards the tube is likely to block due to twisting, making it difficult or impossible to inflate the device.




Step 2:

Insertion

Insert the device using a sterile lubricating cream or gel. The process is similar to inserting a soft vacuum (ventouse) cup.

Make sure that the dome/balloon surface of the device is in contact with the fetal head and the base plate in contact with the pelvic floor.

The device will not inflate or function effectively if placed incorrectly.



Step 3:

Device Position

Once inserted the device should be pushed posteriorly until it is touching the coccyx.

The position is similar to the insertion of a vacuum (ventouse) cup for an occipito posterior position.



Step 4:


Inflation

Patient's legs must be placed flat before inflation is carried out using sterile saline with the 60cc syringe provided.

If the legs are not placed flat before inflating, the device can be expelled or could move during inflation and fail to produce the desired elevation.

A total of 180cc of saline is required to produce the desired elevation (3 syringes of fluid). Close the tap after filling to stop the fluid from leaking.

Inflation volume should not exceed 180cc.



Step 5:

Caesarean Section


Once the inflation is complete, the Caesarean Section is performed using the standard technique.

Step 6:

Device Removal

After delivery of the baby, the two-way tap is opened to release the fluid.

The device is removed by the assistant at the end of procedure by pulling on the tubing or hooking a finger on the plate and pulling the device out of the vagina. If the two-way tap fails, the tube can be cut to release the fluid for removal.



Storage

Store above 5°C and below 30°C. Do not use if package is damaged.

Device Disposal

The device should be discarded according to the hospital regulations.

Clinical studies of fetal pillow use

Randomized control trial of elevation of fetal head with a fetal pillow during caesarean delivery at full cervical dilation¹

This prospective randomized controlled trial was carried out in two teaching hospitals in India and compared the use of Fetal Pillow with other methods of delivery in a second stage Caesarean Section (CS). A total of 240 patients who required a CS in second stage of labor were enrolled into the study. Thirteen patients were excluded from the study, due to lack of informed consent (n=4), previous caesarean (n=2), breech presentation (n=2) and suspected chorioamnionitis (n=5).

Primary Outcome Measure

- Major uterine incision extensions (Grade II and III)

Grade I extensions were defined as those that did not increase operating time and blood loss.

Grade II extensions were defined as those that increased operating time and blood loss.

Grade III extensions were defined as those that involved uterine blood vessels, cervix, vagina or urinary tract.

Secondary Outcome Measures

- Total time taken for CS
- Incision to delivery interval
- Difficulty with delivery of fetal head
- Duration of hospital stay
- Blood loss >1000cc
- Need for blood transfusion
- 5 minute APGAR <3
- NICU stay >24 hours
- Neonatal sepsis
- Neonatal death

Inclusion Criteria

- Ability to give informed consent
- CS at full dilation
- CS after failed instrumental delivery

Exclusion Criteria

- Presence of active genital infection
- Chorioamnionitis
- Breech presentation
- Previous Caesarean Section
- Pregnancy less than 36 weeks
- Inability to give informed consent

Study Methodology

All patients were informed about the trial when admitted to the labor ward. Patients who were able to give informed consent if requiring a CS at full dilation were included in the study. Participants were randomized 1:1 into two parallel groups, the Fetal Pillow group (FP group) and the non-Fetal Pillow group (NFP group).

The delivery methods used in the NFP group were:

- Hand push method
- Reverse breech method
- Abdominal delivery method

CS was carried out using the standard technique and Fetal Pillow was inserted and inflated prior to performing the procedure.

Results

The two groups were similar in terms of their baseline characteristics (Table 1).

Major extensions of uterine incisions were less common in the FP group (Table 2). Incision to delivery interval, total time for CS, need for blood transfusions and length of hospital stay were all significantly lower in the FP group. The intra-operative blood loss >1000mls was more common in the NFP group (Table 2). With regards to fetal outcomes, newborns in the FP group were less likely to have a 5minute Apgar of 3 or less, be admitted to the NICU, or stay in NICU for more than 24 hours than were the newborns in the NFP group (Table 3).