## TYMP

Tympanograph



## otopront



A valuable middle ear diagnosis by means of tympanometry and stapes reflex measurement is standard procedure of an up-to-date ENT office. Otopront's TYMP meets all technical requirements of an up-to-date middle ear analyzer: It works fast, reliable, and is remarkably easy to handle and care in the daily routine use.

A complete automatic measurement procedure can be selected from an easy-to-clean touchscreen which together with a brilliant LCD-display is part of a freely movable frontpanel and thus the device can be integrated in the best possible way into the doctor's ENT unit. The fully automatic test procedure ist started by fitting the probe with the best sealing ear tip into the outer ear canal. Sufficient sealing of the ear tip is indicated by a 2-colour LCD in the probe.

Test procedure is commonly started with tympanometry to get a reliable and objective result on the middle ear pressure and can with the otopront TYMP easily be performed also in restless patients due to its fast operating mode in less than 2 seconds. Stapes reflexes can be recorded ipsi- as well as contra-lateral with different frequencies and loudness levels. In this procedure a pathologic middle ear pressure is automatically compensated by the TYMP and thus provides optimal test conditions. Test results can easily be recorded by the integrated fast printing low-noise thermoprinter and/or transferred via a standard interface (GDT) to the digital patient file.

## **Technical Data**

System Probetone Pressure range

Pressure delta Measure time Reflex tones Sequence of tones Reflex recognition Reflex level Reflexes Pressure calibration State indication:

Change side Interface Air pump: Dimensions Weight Power supply CE Mark:

Accessories:

Microprozessor controlled 226 Hz, 85 dB SPL +200 to -400 daPa Pressure delta: 300 daPa / sec 2 sec. for compliance 500, 1000, 2000 and 4000 Hz at 85, 95 and 105 dB SPL Automatic or manual Automatic Via software adjustable Ipsi and contralateral Automatic 1. LED 3-colors in probe 2. detailed display - symbols Automatic or manual USB 1.1 Very quiet piston pump (W x D x H) 400 x 400 x 140 mm 10 kg 110-220 V~, 50/60 Hz, 60 VA accord. to MDD 93/42/EEC

Tympprobe 30 ear-plugs in 5 different sizes Headset for contralateral test DT 252 Probe -holder

All designs and specifications subject to change without notice.



The swiveling display allows you to handle the TYMP from different angles of view



The printing device can be operated easilily

Happersberger otopront GmbH Langgasse 90 D-65329 Hohenstein Tel. +49 (0) 61 20 - 92 17 0 Fax +49 (0) 61 20 - 92 17 60 www.otopront.com E-Mail:info@otopront.com