Laser Systems



MCO 25plus and MCO 50plus

ALL ESSENTIAL PARAMETERS
CONSTANTLY IN FOCUS



"All essential parameters constantly in focus"

The CO_2 laser is the classic medical laser. It has been used in numerous surgical fields with great success for many years. With the new CO_2 lasers MCO 25plus and MCO 50plus, KLS Martin offers state-of-the-art, cutting-edge laser technology. Whether you consider its uncomplicated operating philosophy, its excellent optical properties or the comprehensive, innovative range of accessories – the MCO25/50plus gives you lots of reasons to decide in favor of a KLS Martin CO_2 laser.

14 good reasons that make a truly convincing case:

- · Novel, futuristic design
- · Quick Tip panel for easy and time-saving parameter selection
- · Easy-moving beam delivery system with excellent beam characteristics
- Safe locking of the articulated (mirror-joint) arm in two different positions
- · Versatile, user-oriented operating modes such as ultra-short single pulses, superpulse mode and cycle mode
- · Reliable beam guidance with minimal spot sizes and high power densities
- Microprocessor control
- Five program storage locations for saving user-defined parameter sets
- Integrated, programmable "Soft Scan plus R" scanner for multi-disciplinary use in conjunction with different handpieces and micromanipulators
- Special LineScan function for robotic vocal cord surgery
- Micro scan patterns for laser-assisted stapedial surgery
- · Micromanipulator Micro Point 2 with minimal spot size of 0.11 mm and highest power densities
- · Possibility to select and turn scan patterns directly via the Micro Point micromanipulator
- · Comprehensive range of accessories for numerous medical fields (ENT, gynecology, dermatology, etc.)



"Take your seat in the cockpit"

Straightforward operating interface, uncomplicated operating philosophy – properties that any user appreciates. The Quick-Tip panel guarantees easy and fast parameter selection. Only the parameters relevant to the operating mode used are automatically displayed at a time – the essentials in focus!



Five program storage locations, enabling you to save your own sets of parameters






"As little as possible, as much as necessary"

The digital scanner system SoftScan plus R comes as an integrated part of the laser unit. It offers controlled, reproducible tissue ablation for use in numerous surgical fields – gentle, tissue-sparing laser application for the benefit of the patient.



for controlled, superficial tissue vaporization in general and aesthetic surgery. Different handpieces and micromanipulators are available for easy handling. The tissue ablation rate can be controlled via different scanning speeds.



Tissue ablation effect



"Unbeatable Variability"

The Micro Point 2 and Mini Point micromanipulators enable you to use the system with all commonly available operating microscopes. Together with the CO_2 laser systems, they constitute the "power pack" for microsurgery – offering extremely high power densities with minimal spot sizes.



76-400-00-04

"Mini Point" micromanipulator for gynecology, proctology Adapters available for connecting operating microscopes and colposcopes from different manufacturers (see page 11)



76-402-00-04

"Micro Point 2" micromanipulator for ENT, neurosurgery (Separate brochure 90-853-02 available on request!) Adapters available for connecting operating microscopes and colposcopes from different manufacturers (see page 11)


"Precision at the highest level"

MCO 25plus and MCO 50plus with scanner system SoftScan plus R and micromanipulator Micro Point 2 R – the highly sophisticated laser system for robotic microsurgery in laryngeal surgery.



76-412-00-04 Micromanipulator "Micro Point 2 R" including Micro Switch for selecting and rotating scan patterns

Cordectomy with LineScan function



Ordering Data

Ordering Data	
KLS Martin CO, Laser MCO 25plus	
consisting of basic unit (25 W),	
articulated (mirror-joint) arm,	
127-mm standard focusing handpiece,	
power cord, foot switch	77-025-01-04
with pilote module ENT	77-025-01-21
KLS Martin CO₂ Laser MCO 25plus	
with integrated scanner	
(standard accessories same as 77-025-01)	77-025-11-04
with pilote module ENT	77-025-11-21
KLS Martin CO₂ Laser MCO 50plus	
consisting of basic unit (50 W),	
articulated (mirror-joint) arm,	
127-mm standard focusing handpiece,	
power cord, foot switch	77-050-01-04
with pilote module ENT	77-050-01-21
KLS Martin CO ₂ Laser MCO 50plus	
with integrated scanner	
(standard accessories same as 77-050-01)	77-050-11-04
with pilote module ENT	77-050-11-21
General accessories	
Focusing handpiece, 50 mm	76 100 05 04
(complete, incl. detachable tip)	76-100-05-04
Tip, detachable,	76 100 06 04
for 50-mm focusing handpiece	76-100-06-04
Focusing handpiece, 127 mm	76-100-10-04
(complete, incl. detachable tip)	76-100-10-04
Tip, detachable,	76-100-11-04
for 127-mm focusing handpiece Focusing handpiece, 200 mm	76-100-11-04
(complete, incl. detachable tip)	76-100-15-04
Tip, detachable,	70-100-13-04
for 200-mm focusing handpiece	76-100-16-04
Backstop for 127-mm	70 100 10 04
focusing handpiece	76-100-20-04
Backstop for 200-mm	
focusing handpiece	76-100-25-04
90° angled working part	
for 127-mm focusing handpiece	76-200-05-04
120° angled working part	
for 127-mm focusing handpiece	76-200-10-04
Adapter for angled working parts	
for 200-mm focusing handpiece	76-200-15-04
CO ₂ laser protective goggles ¹	76-100-50-04
CO ₂ laser protective goggles	
for spectacle wearers ²	76-100-51-04
Combined CO ₂ and Nd:YAG	
laser protective goggles 3	79-100-50-04
Combined CO ₂ and Nd:YAG laser protective	
goggles for spectacle wearers 4	79-100-51-04
Pad for laser protection goggles	
(for item no. 76-100-51-04/79-100-51-04)	79-101-51-04
Strap for laser protection glasses	
(for item no. 76-100-51-04/79-100-51-04)	79-101-50-04
All laser units are available with or without scal	


Ordering Data	
Micromanipulators	
"Mini Point" (gynecology, proctology)	76-400-00-04
"Micro Point 2" (ENT, neurosurgery) 5	76-402-00-04
"Micro Point 2 R" (incl. Micro Switch)	76 410 00 04
(ENT, neurosurgery) ⁵	76-412-00-04
Adapters for Mini Point	
Adapter, Zeiss system	76-400-10-04
Adapter, Möller-Wedel system	76-400-12-04
Adapter, Leica system	76-400-13-04
Adapter, Zeiss OPMI FC 1	76-400-14-04
Adapter, Möller-Wedel VM 500	76-400-15-04
Adapter, Leica M 841	76-400-16-04
Adapter, Olympus OCS-500	76-400-17-04
Adapter, Zeiss 150 FC	76-400-18-04
(other adapters available on request)	
Adapters for Micro Point 2 and Micro Point 2 R	
Adapter, Leica M 650/651/655 systems	
(other adapters available on request)	76-401-01-04
Adapter, Leica M 680	76-401-02-04
Adapter, Möller-Wedel VM 900	76-401-02-04
Adapter, Moher-Weder VM 900 Adapter, Leica M 841, complete	76-401-16-04
Adapter, Olympus OCS-500	76-401-17-04
Adapter, Opmi Drape, complete	76-401-17-04
Adapter, Opmi 6, complete	76-401-19-04
Adapter, Leica M 520, complete	76-401-19-04
Adapter, Leica W 320, complete	70-401-20-04
Adapters for Mini Point, Micro Point 2 and Micro I	Point 2 R
Adapter, Kaps system 42/52/62	76-401-11-04
Micro Switch for micromanipulator Micro Point 2 in connection with LineScan Scanner	76-401-30-04
"SoftScan plus R" scanner retrofit kit	
for MCO 25plus	76-500-01-04
"SoftScan plus R" scanner retrofit kit	76 500 00 04
for MCO 50plus	76-500-02-04
The retrofit kits (complete) include:	76 500 21 04
- scanner head, complete	- 76-500-31-04
- scanner connecting cable	- 76-500-40-04
- scanner handpiece, 127 mm	76 500 50 04
(complete, incl. tip)	- 76-500-50-04
- scanner handpiece, 200 mm (complete, incl. tip)	- 76-500-60-04
- incl. all technical components	- 76-300-60-04
Adapter for connecting "SoftScan plus R"	
to "Mini Point"	76-500-41-04
Adapter for connecting "SoftScan plus R"	76-300-41-04
	76-500-42-04
to "Micro Point"	
Tip for 127-mm handpiece	76-500-51-04 76-500-61-04
Tip for 200-mm handpiece	76-300-61-04
Special accessories	
Special optics set for nasal operations,	76-600-00-04
consisting of:	
- optics	- 76-600-01-04
- 0-degree attachment for optics	- 76-600-02-04
- 90-degree attachment for optics	- 76-600-03-04
- 120-degree attachment for optics	- 76-600-04-04
- brush for attachments	- 76-600-05-04
 D 9000–10600 L5 RH DIN CE acc. to 89/686/EEC (PF D 9000–11000 L4 RH DIN CE acc. to 89/686/EEC (PF 	
D 9000–11000 L4 RH DIN CE acc. to 89/886/EEC (PPE) /	L/
DI 10600 nm L4 RH DIN CE	
D 1060 L7 RH DIN CE acc. to 89/686/EEC (PPE) / DI 10600 nm L3 RH DIN CE	
⁵ ZEISS adapter inclusive	

Handpieces and Attachments




* included in Laser package

** included in Laser with Scanner package



76-200-15-04 Adapter for angled working parts for 200-mm focusing handpiece 76-200-05-04 90° angled working part, for 127-mm focusing handpiece and 200-mm focusing handpiece (in conjunction with adapter 76-200-15-04) 76-200-10-04
120° angled working part,
for 127-mm focusing handpiece
and 200-mm focusing handpiece
(in conjunction with adapter 76-200-15-04)

76-600-00-04 Special optics set for nasal operations, consisting of:



76-600-05-04 Brush for ENT attachments

Technical Data - MCO 25plus				
Type of laser	MCO 25plus			
Wavelength	10.6 μm, infrared			
Laser power (applied to tissue)	2-25 W			
CO ₂ laser tube	sealed off			
	DC-excited			
Operating modes:				
continuous-wave (CW)	2–25 W			
superpulse	0.3 ms; 11 W, average, max. power			
single pulse	5 ms – 10 s			
omgre pane	output power: 10–25 W			
pulse train	minimum pulse length: 5 ms			
puise train	average output power: < 0.1 – 25 W			
avaliably repeated pulses				
cyclically repeated pulses	cycle length: - permanent			
	- 10 ms – 10 s			
Dilat Iaaan liaht	Diede lees entieve die die tele			
Pilot laser light	Diode laser, continuously adjustable			
	635 nm, ≤ 5 mW			
	light-red			
Laser beam delivery	Articulated spring arm with 7 joints/mirrors,			
	arm length = 1,300 mm,			
	handpiece exchangeable			
Focus diameter	- 0.15 mm with a focal distance			
	of 50 mm (spot size)			
	- 0.20 mm with the standard focal distance			
	of 127 mm			
	- 0.30 mm with a focal distance			
	of 200 mm			
Operator/control panel	Bluemode LCD (liquid crystal display)			
	with illuminated buttons featuring symbols,			
	keypad sealed, with integrated scanner			
	functions			
Program memory	5 freely assignable memory locations			
	(can be used for scanner as well)			
Cooling	internal			
Power requirements	230 V / 4 A			
·	50 / 60 Hz			
Protection class	1			
Type (of applied part)	В			
Laser class	4			
Dimensions	290 x 450 x 1,950 mm			
	(height with articulated arm folded down:			
	1,180 mm)			
Weight	57 kg			
HOISIIL				



Technical Data - MCO 50plus

Type of laser	MCO 50plus
Wavelength	10.6 μm, infrared
Laser power (applied to tissue)	2-50 W
CO ₂ laser tube	sealed off DC-excited
Operating modes:	
continuous-wave (CW)	2-50 W
superpulse	0.3 ms; 16 W, average, max. power
single pulse	10 ms - 10 s output power: 10-50 W
pulse train	minimum pulse length: 8 ms average output power: < 0.1 - 50 W
cyclically repeated pulses	cycle length: - permanent - 10 ms - 10 s
Pilot laser light	Diode laser, continuously adjustable 635 nm, \leq 5 mW light-red
Laser beam delivery	Articulated spring arm with 7 joints/mirrors arm length = 1,300 mm, handpiece exchangeable
Focus diameter	 - 0.15 mm with a focal distance of 50 mm (spot size) - 0.20 mm with the standard focal distance of 127 mm - 0.30 mm with a focal distance of 200 mm
Operator/control panel	Bluemode LCD (liquid crystal display) with illuminated buttons featuring symbols keypad sealed, with integrated scanner functions
Program memory	5 freely assignable memory locations (can be used for scanner as well)
Cooling	internal
Power requirements	230 V / 4 A 50 / 60 Hz
Protection class	I
Type (of applied part)	В
Laser class	4
Dimensions	290 x 450 x 2,170 mm (height with articulated arm folded down: 1,630 mm)
Weight	62 kg





Technical Data - Optional Accessories

"Mini Point" micromanipulator					
Zoom range continuously adjustable between 200 and 400 mm					
Focal length	200 mm	250 mm	300 mm	350 mm	400 mm
Spot diameter	0.36 mm	0.47 mm	0.56 mm	0.65 mm	0.76 mm
Power density (with 25 W)	24.5 kW/cm ²	14.5 kW/cm ²	10.2 kW/cm ²	7.5 kW/cm ²	5.5 kW/cm ²

"Micro Point 2" / "Micro Point 2 R" micromanipulator							
Zoom range continuously adjustable between 225 and 500 mm							
Focal length	225 mm	250 mm	300 mm	350 mm	400 mm	450 mm	500 mm
Spot diameter	0.110 mm	0.120 mm	0.130 mm	0.135 mm	0.140 mm	0.145 mm	0.150 mm
Power density (with 25 W)	263 kW/cm ²	221 kW/cm ²	188 kW/cm ²	175 kW/cm ²	163 kW/cm ²	151 kW/cm²	141 kW/cm ²

"SoftScan plus R" scanner

 $\ensuremath{\mathsf{6}}$ different scan patterns in different sizes and positions for tissue vaporization

3 different line patterns with variable length for robotic tissue resection

Electronically rotatable micro scan patterns ranging from 0.4 mm to 2.0 mm, especially for otology

Rotation via micromanipulator Micro Point R

Operating/control functions integrated into the laser unit's display screen

Irradiation time per treatment spot (depending on selected parameters) 0.15--1.2~ms

 $\label{thm:constraint} \mbox{Scan pattern representation with pilot laser during treatment pauses as well}$

Energy densities with 127-mm handpiece: MCO 25plus: max. 24 J/cm^2


MCO 50plus: max. 48 J/cm²

Innovative Laser Systems from KLS Martin

For further information on our laser units, please contact:

Gebrüder Martin GmbH & Co. KG A company of the KLS Martin Group Ludwigstaler Str. 132 · D-78532 Tuttlingen Postfach 60 · D-78501 Tuttlingen Tel. +49 7461 706-0 · Fax +49 7461 706-193 info@klsmartin.com · www.klsmartin.com







KLS Martin MCO 25plus

KLS Martin Group

Karl Leibinger GmbH & Co. KG

78570 Mühlheim Germany Tel. +49 7463 838-0 info@klsmartin.com

Rudolf Buck GmbH

78570 Mühlheim Germany Tel. +49 7463 99516-30 info@klsmartin.com

Martin Nederland/Marned B.V.

1270 AG Huizen The Netherlands Tel. +31 35 523 45 38 nederland@klsmartin.com

KLS Martin L.P.

Jacksonville, FI 32246 USA Office phone +1 904 641 77 46 usa@klsmartin.com Stuckenbrock Medizintechnik GmbH

78532 Tuttlingen Germany Tel. +49 7461 161114 verwaltung@stuckenbrock.de

KLS Martin France SARL

68000 Colmar France Tel. +33 3 89 21 66 01 france@klsmartin.com

Nippon Martin K.K.

Osaka 541-0046 Japan Tel. +81 6 62 28 90 75 nippon@klsmartin.com

Orthosurgical Implants Inc.

Miami, FI 33186 USA Office phone +1 877 969 45 45 sales@orthosurgical.com KLS Martin GmbH + Co. KG

79224 Umkirch Germany Tel. +49 7665 98 02-0 info@klsmartin.com

Martin Italia S.r.l.

20059 Vimercate (MB) Italy Tel. +39 039 605 67 31 italia@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Representative Office 121471 Moscow Russia Tel. +7 (499) 792-76-19 russia@klsmartin.com

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group Ludwigstaler Str. 132 · D-78532 Tuttlingen Postfach 60 · D-78501 Tuttlingen Tel. +49 7461 706-0 · Fax +49 7461 706-193 info@klsmartin.com · www.klsmartin.com

