



VELNEZ[®]

The Optimal Healing Solution

VELNEZ®

The Optimal Healing Solution



CE
2460

VelNez® is a fragmentable composite that fragments within a few days on application intended as a nasal dressing. The pain associated with the traditional nasal pack removal is completely eliminated with VelNez® since such a procedure is not necessary as VelNez® fragments uniformly. It will also promote healing and reduce the clotting time.

Intended Use

The fragmentable VelNez® is intended to be used in the nasal/ sinus cavity in surgery or trauma patients. It acts as a space occupying dressing, preventing adhesion by separating the compromised mucosal surfaces. It also helps in healing and achieving haemostasis.

Instructions for Use

Step 1: Removal from Package - Open the blister pack by removing the paper film and take out VelNez®.

Step 2: Pre-Insertion – if required cut the VelNez® using sterile apparatus.

Step 3: Insertion - Using sterile nasal dressing forceps place the VelNez® inside the nasal cavity at the site of application and irrigate using saline solution if required.

Step 4: Confirm - Confirm that VelNez® is placed at the intended position. VelNez® will fill itself by absorbing the surrounding fluids or blood.

*Note: The VelNez® dressing may be irrigated with saline water as and when required.

Constituents

Gelatin*(BSE/TSE agentfree) Chitosan, Polyvinyl Alcohol, Psyllium Husk.

Warning & Precautions

- Do not re-sterilize.
- Do not re-use.
- Use the device prior to the "Use by (📅)" specified on package.
- To be used by Registered Medical/ Practitioner or a hospital.
- Do not use if package is damaged.
- Discontinue the use if the application area shows signs of infection, irritation, maceration (whitening of surrounding skin) or itching & consult a healthcare professional.
- Should not be used for the patients with Coagulation disorders.

Risk of Reuse

VelNez® is supplied sterile for single use only. Do not re-sterilize as this may change product characteristics and lead to failure of Intended use of the product.

Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and /or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease from one patient to another and loss of uptake of device and reduced the efficacy is prominent. No other consequences will happen.

Damaged, unused, opened packages of VelNez® should be discarded as per disposal instructions, sterility of the product is not guaranteed when the packaging is damaged or unused portion of the device is left in the pack.

Indications

As nasal pack for hemorrhage control or/and prevention of adhesion after nasal surgery, Septoplasty with and without Functional endoscopic sinus surgery (FESS) and turbinoplasty.

Contraindications

No known contraindications and adverse reactions reported.

Note: Do not use VelNez® if you are allergic (Hypersensitive) to its constituents.

Disposal of Device

Discard as per hospital biomedical waste management system prevailing in the premises meeting local or national regulatory laws.

Patient Population

Patients with age group from 18-60 years of age.

Intended User

Surgeon, ENT specialist

How Supplied

VelNez® is supplied as a sterile, single-use Nasal pack.

Storage Conditions

VelNez® should be stored at temperature between 4°C to 30°C. Keep away from direct sunlight. Store in Dry place.

Shelf Life

3 years from the date of Mfg.

Method of Sterilization

Gamma irradiation

Available sizes

Shape	Size
Rectangle	4 cm x 2 cm
Rectangle	8 cm x 2 cm

PATENTED

For more information write to Marketing Department



Datt Mediproducts Pvt. Limited

(An ISO 13485 Certified Company)

Admin Office: 56, Community Centre, East of Kailash, New Delhi-110065

Works: 52-54 & 63-64, Roz-ka-Meo Industrial Area, Nuh, Distt. Mewat-122 103, Haryana, India

E-mail: advancedevice@dattmedi.com | Website: www.dattmedi.com

Customer Care Toll Free No.: 1800-1030-707

ECREP

mdi Europa GmbH Langenhagener Str. 71

D-30855Hannover-Langenhagen Germany Phone: +49-511-3908 9530

Internet: www.mdi-europa.com Trade Register Hanover HRB 58423