


Inflatable tampanade designed to create a fast, single treatment solution with a low-profile and self-lubricating properties




RAPID RHINO° Epistaxis Products

Offering a full range of epsitaxis solutions




Features

- Easy to insert and remove
- Quickly conforms to nasal anatomy
- Provides gentle and even compression to areas of epistaxis
- Designed to prevent adherence to tissue or blood clots upon removal
- Available in a wide range of sizes



RAPID RHINO Nasal Pac (Uninflated 5.5cm anterior shown)



RAPID RHINO with Carboxymethylcellulose (CMC)

Comparison


A comparison of RAPID RHINO° epistaxis products to Rhino Rocket² for the treatment of epistaxis in the emergency department¹: RAPID RHINO epistaxis products provide greater patient comfort at insertion and at removal.




RAPID RHINO product usage directions




1 Soak in sterile water for a FULL 30 seconds.




2 Insert along superior aspect of the hard palate until the blue indicator is past the nares.




3 Using a 20ml syringe, inflate the Rapid Rhino device with AIR only. Monitor the pilot cuff for direct tactile feedback; Stop inflation when the pilot cuff becomes rounded and feels firm when squeezed.



4 Inflate the cuff to provide a gentle, lowpressure tamponade delivering the CMC fabric directly to the bleed site.



5 Reassess after 15-20 minutes; reinflate to ensure proper pressure (if necessary) and tape to patient's cheek away from the upper lip.



6 Removal should occur 24-72 hours after treatment.

Unilateral models



R450 nasal tamponade – 4.5cm inflatable epistaxis tamponade device



RR 530 non-inflatable tamponade RAPID-PAC° expistaxis tamponade device is 5.5cm offers a noninflatable option for epistaxis



RR550 inflatable tamponade – 5.5cm inflatable epistaxis tamponade device



RR551 inflatable tamponade with airway – 5.5cm inflatable epistaxis tamponade device that contains an airway



RR750 inflatable tamponade – 7.5cm inflatable epistaxis tamponade device



RR751 inflatable tamponade with airway – 7.5cm inflatable epistaxis tamponade device that contains an airway



RR900 – 9cm inflatable tamponade device with two balloons to control intranasal hemorrhage originating from the sphenopalatine artery

Bilateral models



RR552 5.5cm bilateral epistaxis device; RR555 5.5cm bilateral epistaxis device with an airway



RR752 7.5cm bilateral epistaxis device; RR755 7.5cm bilateral epistaxis device with an airway

References

- 1 Singer AJ, Blanda M, Cronin K, LoGiudice-Khwaja M, Gulla J, Bradshaw J, Katz. A Comparison of nasal tampons for the treatment of epistaxis in the emergency department: a randomized controlled trial. Ann Emerg Med. 2005; 45:134–139
- 2 Rhino Rocket is a registered trademark of Shippert Medical Technologies Corp. and is used for reference purposes only.

Ordering information

RAPID RHINO° inflatable tamponade devices

Reference #	Description
RR450	4.5CM INFLATABLE UNILATERAL WITHOUT AIRWAY
RR530	5.5CM NON-INFLATABLE RAPID-PAC°
RR550	5.5CM INFLATABLE UNILATERAL WITHOUT AIRWAY
RR551	5.5CM INFLATABLE UNILATERAL WITH AIRWAY
RR552	5.5CM INFLATABLE BILATERAL WITHOUT AIRWAY
RR555	5.5CM BILATERAL EPISTAXIS DEVICE WITH AIRWAY
RR750	7.5CM INFLATABLE UNILATERAL WITHOUT AIRWAY
RR751	7.5CM INFLATABLE UNILATERAL WITH AIRWAY
RR752	7.5CM INFLATABLE BILATERAL WITHOUT AIRWAY
RR755	7.5CM INFLATABLE BILATERAL WITH AIRWAY
RR900	9CM INFLATABLE DUAL BALLOON

Smith & Nephew, Inc.

7000 West William Cannon Drive

Austin, TX 78735 USA www.smith-nephew.com

Information: 1-800-343-5717 Orders/Inquiries: 1-888-994-2782