

## MVR® Mount

The MVR® Mount is especially designed for MediCap® MVR® Medical Video Recorders. The recorder can be assembled directly on your cart to offer maximum flexibility: vertically, horizontally, left or right.



### Flexibility

Space-saving design allows for maximum cart usage.

### Secure and Easy Assembly

Mounting kit comes with everything for fast, secure mounting.


### Cable Management

Benefit from extra space for the cable arrangement.

Medical Imaging Made Easy™






## Examples of mounting



### Direct assembly to a column

Use your existing space by assembling your MVR® recorder vertically inside or outside a column.




### Bracket for horizontal assembly

The additional bracket offers even more flexibility. Assemble your MVR® Recorder in any suitable position.

<sup>1</sup> More examples on [www.medicapture.com](http://www.medicapture.com)

## MVR® Mount Specs<sup>2</sup>



### Easy Assembly

Use the relocatable nut stones of your video cart or any other fixation to assemble the MVR Mount like your cart accessories.

**Order number** MNT0001MVR

**Material** Aluminium  
Powder coated, RAL 9016  
(white)

**Mechanical** Main mount:  
Size (l x w x h): 215 x 168 x 48mm  
/ 8.46 x 6.61 x 1.89"  
Retaining bracket:  
Size (l x w x h): 165,9 x 110 x 77mm  
/ 6.53 x 4.33 x 3.03"  
Weight: 0,7kg / 1.54lb

**Certifications** RoHS

**Kit Contents** MediCap® MVR® Mount  
(pre-assembled),  
Mount retaining bracket incl. screws,  
M4 nut stones for 8mm rails, M4  
screws, M6 screws, Screwdriver, Alan  
key, Installation Guide.

<sup>2</sup> Subject to change.



### MediCapture, Inc. – North America

2250 Hickory Road, Suite 200  
Plymouth Meeting, PA 19462, USA

+1-888-922-7887

[www.medicapture.com](http://www.medicapture.com)  
[sales@medicapture.com](mailto:sales@medicapture.com)

### MediCapture Europe B.V.

Kraanpoort 10  
6041 EG Roermond  
The Netherlands

+31 475 214 305

[www.medicapture.com](http://www.medicapture.com)  
[sales.europe@medicapture.com](mailto:sales.europe@medicapture.com)

### MediCapture APAC Inc.

3FL-6, No. 432, Section 1, Keelung Road,  
Xinyi District 110, Taipei City,  
Taiwan, Republic Of China

+886- 2- 2345- 3668

[www.medicapture.com](http://www.medicapture.com)  
[sales.asiapacific@medicapture.com](mailto:sales.asiapacific@medicapture.com)

© 2017 MediCapture, Inc.  
MVRMOUNT-DS01-EN-170131-MC

Medical Imaging Made Easy™