

A New Convenient System



Disposable Proctorectoscope Tubes

Packed sterile – for immediate use

Disposable Proctorectoscope Tubes

A New Convenient System

Richard Wolf Disposable Proctorectoscope Tubes have a long track record of providing a clear view in proctological diagnostics for many years.

This step has now become even more efficient. Sterile packing at the factory means that the proctorectoscope tubes are now ready for use immediately. This eliminates the process of disinfection prior to application.



■ Improved handling

Each tube is packed singly and sterile. Disinfection prior to application is eliminated!

■ Cost-effective solution

You benefit from the ideal combination of disposable tube and reusable handle.

■ Everything in view

An easily readable scale makes documentation and checks easier.

■ Everything under control

The reusable handle permits easy application with optimum illumination.

Tubes and accessories

	Rectoscopy tube , sterile, Ø 21 mm, working length 255 mm, pack of 104844255
	Proctoscopy tube , sterile, Ø 21 mm, working length 120 mm, pack of 104844120
	Handle with annular cold light illumination8838.90
	Hand bulb, sterile (pack of 5) Hand bulb with filter for disposable use Ø 49 mm, OL 567 mm412700
	Locking cap with observation window8838.60