

# AS300

innovative, reliable, adaptive



**Soluscope**

by  Laboratoires  
**ANIOS**  
Le professionnel de la désinfection

# AS300



an

# Efficient

concept



AS300 pass through

*A high-level security for your endoscopes*

#### INNOVATIVE AIR FLOW SYSTEM

The patented system of laminar air flow filtered with an H14 high capacity filter protects your endoscopes from any contamination.

The overpressure constantly maintained in the cabinet prevents any entry of airborne micro-organisms. The high rate of air renewal maintains the microbiological quality of your endoscopes.

#### SERVICE WORKFLOW OPTIMIZATION

Your endoscopes are ready-to-use at any time. The automatic carousel provides direct access to the endoscopes. It reduces handling manipulations and eliminates identification faults.

Each endoscope is plugged to the cabinet with a dedicated connector that can be easily reprocessed in *Soluscope AERs* to increase your productivity in total safety.

CE MARK / EN 16442\* compliant

# a Scalable device



AS300 one-door (single/double)

## *Designed for evolving needs*

### **FLEXIBLE AND MODULAR**

The AS300 allows vertical storage for 10 endoscopes validated up to 7 days. An optional compartment can be simply integrated to extend the AS300 capacity to 20 endoscopes. To reinforce hygiene of your endoscopes and optimize the workflow, the AS300 is available in a pass through version in order to physically separate the white room from the grey room. For enhanced safety of your endoscopes, accessories are available to prevent endoscope tangling and to hold ultrasound connectors. Individual compartments for endoscope valves are provided to secure their storage.

### **A COMPLETE TRACEABILITY**

The embedded traceability system includes barcode or RFID reader and printer. Traceability tickets are printed out and securely stored on the AS300 memory card. The Ethernet port enables connection to the healthcare facility network and is compliant with the *IT Soluscope* data management system (optional).

storage cabinet for endoscopes

a  
**Clever**  
approach



*A solid and effortless device*

**DESIGNED FOR MEDICAL USE**

The AS300 offers a user friendly touch screen interface with a colour code as a decision-making aid.

A spherical shape with smooth surfaces and a wide door opening facilitate cleaning.

**A ROBUST MACHINE**

The AS300 is a cast solid stainless steel cabinet and offers a high level of finishing.

The quality of cabinet components guarantees a minimal maintenance facilitated by an easy and secure access to technical compartments.

# A complete solution to optimize your organization

## TRANSPORT

*Soluscope Trolley* is designed to dispatch endoscopes from treatment to examination room. The specific tray equipped with a transparent lid safely stores the endoscope during its transportation. *Soluscope Trolley* facilitates your daily routine in a secure manner.

## STORAGE

*Soluscope Rack* can store up to 10 trays in a reduced space.

## SERVICES

After-sales service, package and training are designed to get the best of your device.

## TECHNICAL SPECIFICATIONS

Validated storage time	Up to 7 days		
Electrical supply	220-240V/50Hz		
Air supply	Connection to medical air		
Air treatment exterior endoscope surface	Prefilter and HEPA filter H14		
Air treatment in the channels	Medical air		
Traceability	Barcode reader - Printer - Removable memory card - RFID Reader (option) - IT Soluscope (option**)		
Connector treatment	Soluscope AER disinfection or sterilisation (autoclave)		
Version	Single	Double	Pass through
Capacité (endoscopes)	10	20	10
Dimensions en cm (L x P x H)	116x73x235	184x73x235	110x80x217
Poids	300 kg	450 kg	286 kg



\* validated by Eurofins Biotech-Germande  
 \*\* coming soon





## SOLUSCOPE COMPLETE SOLUTION

*Soluscope* proposes a complete solution for flexible endoscope reprocessing.

It offers products for every step of the reprocessing cycle.

Each part is certified, easy to use and effective in order to reduce your daily load.



[www.soluscope.com](http://www.soluscope.com)

Soluscope SAS : 100, rue du Fauge - Z.I. Les Paluds - 13400 Aubagne - France

☎ +33 (0)4 91 83 21 22 📠 +33 (0)4 91 83 21 10 ✉ info@soluscope.com

*Soluscope* products are compliant with *European and international norms*, and are *CE marked*. *ISO 9001 version 2008 and ISO 13485 version 2003 certified*, *Soluscope* is continuously committed to satisfying its clients rigorous demands.

distributed by

