




Innovation in Gynecology



LiNA Pencils

*For open monopolar
electrosurgery*



LiNA Surgipen Reusable



With integrated silicone cable. Autoclavable 100 times. Available with standard US-3-pin connector. Available for 2.35mm electrodes and for 4mm hexagonal electrodes. With 3 or 5m cable.



REF 508/512

LiNA Pen, fingerswitch controlled Single Use

With 3m cable and US-3-pin connector. Available for 2.35mm electrodes.



REF 516

Ordering information:

LiNA Pencils

REF:	Product	Packing
508-L-3	LiNA Surgipen - Reusable - Autoclavable - for 2.35mm electrodes. With 3m silicone cable and standard US-3-pin connector. With 70mm blade electrode.	Box of 10, non-sterile
508-L-5	LiNA Surgipen - Reusable - Autoclavable - for 2.35mm electrodes. With 5m silicone cable and standard US-3-pin connector. With 70mm blade electrode.	
512-L-3	LiNA Surgipen - Reusable - Autoclavable - for 4mm hexagonal electrodes. With 3m silicone cable and standard US-3-pin connector.	Box of 10, non-sterile
512-L-5	LiNA Surgipen - Reusable - Autoclavable - for 4mm hexagonal electrodes. With 5m silicone cable and standard US-3-pin connector.	
516	LiNA Pen, fingerswitch controlled – Single use – for 2.35mm electrodes. With 3m cable and standard US-3-pin connector. With 70mm blade electrode.	Box of 25, sterile
516-NS	LiNA Pen, fingerswitch controlled – Single use – for 2.35mm electrodes. With 3m cable and standard US-3-pin connector. With 70mm blade electrode.	Box of 25, non-sterile

Quality Management System
NF EN ISO13485:2016