

# Easy-care through to the OR



**Maintenance-free  
Two-way Stopcocks**

# Maintenance-free Two-way Stopcocks

Easy-care through to the OR

We have now eliminated grease and oil in the maintenance-free two-way stopcocks so as to obtain a maximally reliable result during sterilization.

The new stopcock can also be dismantled and assembled by simply "snapping" into place: This leaves you more time for the important things.



- Easy dismantling and assembly
- No grease or oil necessary – completely maintenance-free!
- Unrestricted machine reprocessing
- Sterilization using all standard procedures and in assembled condition!

## Suction-irrigation tube

with maintenance-free two-way stopcock,  
sheath tube with distal openings, reusable

WL 290 mm ..... **8583732**  
WL 410 mm ..... **8584732**

## Suction-irrigation adapter

with maintenance-free two-way stopcock,  
for multifunctional suction-irrigation system,  
reusable

..... **8583819**

## Two-way stopcock maintenance-free

with luer-lock connector, reusable

Inside diameter 4 mm ..... **88351**  
Inside diameter 3 mm ..... **88350**

