


Rhinology

Expanding PVA Nasal Tampons, Anatomical

- Anatomical form
- Post-operative tamponing after rhinological procedures
- Silicone airway (if required – for easier breathing and swallowing)
- 15 cm safety thread



IMPORTANT!

Period of usage should not exceed 4 days.

Art. No.	Description	Length	Packaging
66140	Nasal tampon with airway	75 mm	20 pieces / box
66150	Nasal tampon	75 mm	20 pieces / box
66160	Nasal tampon, slim	75 mm	20 pieces / box

Expanding PVA Epistaxis Tampons

- Application as epistaxis tampon to arrest bleeding in medial or posterior nasal area, or post-operative usage
- Three sizes available
- 15 cm safety thread



Art. No.	Size	Length	Packaging
66200	Epistaxis tampon large	100 mm	20 pieces / box
66210	Epistaxis tampon medium	55 mm	20 pieces / box
66220	Epistaxis tampon pediatric	35 mm	20 pieces / box

Expanding PVA Sinus Tampons

- Post-operative application after sinus surgical procedures
- Special shape for tamponing of the ethmoid cavity
- 15 cm safety thread




Art. No.	Description	Dimensions	Packaging
66300	Sinus tampon	35 x 9 x 12 mm	20 pairs / box
66310	Sinus tampon slim	35 x 6 x 12 mm	20 pairs / box



Foam Nasal Packing with Latex-free Non-adhesive Cover

- Post-operative application for rhino-/ septoplasty
- Water repellent and non-adhesive
- Safety fixation of thread
- X-ray contrast strip
- Sterile packed in pairs
- 25 cm safety thread


Art. No.	Dimensions	Packaging
35083	70 x 20 x 10 mm	25 pairs / box
35084	50 x 20 x 10 mm	25 pairs / box
35085	80 x 20 x 10 mm	25 pairs / box
35086	70 x 10 x 10 mm	25 pairs / box



Foam Sinus Packing with Latex-free Non-adhesive Cover


- Post-operative application for rhino-/ septoplasty
- Water repellent
- Non-adhesive
- Safety fixation of thread
- X-ray contrast strip
- Sterile packed in pairs
- 25 cm safety thread

Art. No.	Dimensions	Packaging
35087	35 x 10 x 8 mm	25 pairs / box
35088	35 x 10 x 6 mm	25 pairs / box
35089	25 x 10 x 6 mm	25 pairs / box
35091	35 x 10 x 10 mm	25 pairs / box



Foam Nasal Packing with Latex Cover

- Application after rhino-/ septoplasty
- Non-adhesive
- Easy application
- 30 cm safety thread
- X-ray contrast strip
- 16 different sizes available



Art. No.	Dimensions	Packaging
35043	70 x 20 x 10 mm	25 pairs / box
35044	50 x 20 x 10 mm	25 pairs / box
35045	60 x 20 x 10 mm	25 pairs / box
35046	80 x 20 x 10 mm	25 pairs / box
35048	80 x 20 x 20 mm	25 pairs / box
35049	80 x 20 x 15 mm	25 pairs / box
35050	50 x 20 x 15 mm	25 pairs / box
35051	60 x 20 x 15 mm	25 pairs / box

Art. No.	Dimensions	Packaging
35052	70 x 20 x 15 mm	25 pairs / box
35053	50 x 20 x 20 mm	25 pairs / box
35054	60 x 20 x 20 mm	25 pairs / box
35055	70 x 20 x 20 mm	25 pairs / box
35057	50 x 10 x 10 mm	25 pairs / box
35058	60 x 10 x 10 mm	25 pairs / box
35059	70 x 10 x 10 mm	25 pairs / box
35069	80 x 10 x 10 mm	25 pairs / box