







Linear stapler

Unique safety feature!

The color indicator window shows different colors for different firing phases and reduces the possibility of misoperation.



-  *primary status*
-  *semi-closed status*
-  *fully-closed status*
-  *post-firing status*

PRODUCT INFORMATION

<i>Product Code</i>	<i>Description</i>
FLSLF30 3.5mm	Linear stapler with color indicator 30mm, 3.5mm staple height
FLSLF30 4.8mm	Linear stapler with color indicator 30mm, 4.8mm staple height
FLSLF45 3.5mm	Linear stapler with color indicator 45mm, 3.5mm staple height
FLSLF45 4.8mm	Linear stapler with color indicator 45mm, 4.8mm staple height
FLSLF60 3.5mm	Linear stapler with color indicator 60mm, 3.5mm staple height
FLSLF60 4.8mm	Linear stapler with color indicator 60mm, 4.8mm staple height
FLSLF90 3.5mm	Linear stapler with color indicator 90mm, 3.5mm staple height
FLSLF90 4.8mm	Linear stapler with color indicator 90mm, 4.8mm staple height
(FLSL)H3.5A	Reloading Unit for Linear Stapler 30mm, 3.5 mm staple height
(FLSL)H4.8A	Reloading Unit for Linear Stapler 30mm, 4.8 mm staple height
(FLSL)H3.5B	Reloading Unit for Linear Stapler 45mm, 3.5 mm staple height
(FLSL)H4.8B	Reloading Unit for Linear Stapler 45mm, 4.8 mm staple height
(FLSL)H3.5C	Reloading Unit for Linear Stapler 60mm, 3.5 mm staple height
(FLSL)H4.8C	Reloading Unit for Linear Stapler 60mm, 4.8 mm staple height
(FLSL)H3.5E	Reloading Unit for Linear Stapler 90mm, 3.5 mm staple height
(FLSL)H4.8E	Reloading Unit for Linear Stapler 90mm, 4.8 mm staple height

The information presented in this material is intended for product description and promotion only. Panther Healthcare has the relevant power of interpretation. Always refer to the product label and/or instruction for use before using any Panther Healthcare product. For further details and updates please refer to Panther Healthcare website: www.pantherhealthcare.com