## **CLEANASCOPE™ ADVANTAGE**

Transport & Short Term Storage System



# CLEAR IDENTIFICATION OF CLEAN VERSUS CONTAMINATED ENDOSCOPES



### CONFIDENCE IS KNOWING IT'S CLEAN AT A GLANCE

The CLEANASCOPE<sup>™</sup> Advantage Transport and Short-term Storage System protects your endoscopes and other delicate equipment by minimising direct handling and reducing contact with airborne contaminants.

Its patented design creates a protective environment for each disinfected endoscope. Green (clean) and red (contaminated) colour-coded liners immediately identify clean or contaminated equipment.

## OPTIMISED ENDOSCOPE WORKFLOW MAXIMISES INFECTION PREVENTION



Line tray using the sterile base liner, and place the clean endoscope and the folded red liner into the tray.





REPROCESSED ENDOSCOPE

Cover the tray with the sterile Green liner marked "CLEAN" identifying this as a patient ready endoscope.



STORE OR TRANSPORT PATIENT READY ENDOSCOPE



#### THE COMPLETE CIRCLE OF PROTECTION



As the global vanguard in infection prevention, **only Cantel delivers the Complete Circle of Protection**, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimize your success.

DRY & STORAGE Bacteria pose significant risk to endoscopes during transport and storage. Cantel's transport, drying and storage solutions are designed to protect valuable inventory, reduce cross contamination touchpoints, eliminate moisture in the endoscope channels and control humidity. Humidity or moisture in endoscope channels is known to aid bacterial growth.





**CONTAMINATED TRANSPORT** 





MANUAL CLEANING



CLEANASCOPE™ CART ORDERING INFORMATION			
PART NUMBER	PRODUCT DESCRIPTION	DIMENSIONS	
103470	Five tray cart (trays included)	88.14H x 64.00W x 56.13D (cm)	
103466	Six tray cart (trays included)	102.87H x 64.00W x 56.13D (cm)	
103467	Ten tray cart (trays included)	161.29H x 64.00W x 56.13D (cm)	







"Transport the soiled endoscope to the reprocessing area in a closed container that prevents exposure to staff, patients or the environment to potentially infectious organisms"

ASGE Quality Assurance in Endoscopy Committee, 2011 and SGNA/ESGENA Standards of Infection Control in Reprocessing of Flexible Gastrointestinal Endoscopes 2012

CLEANASCOPE™ LINERS & ACCESSORIES ORDERING INFORMATION			
PART NUMBER	PRODUCT DESCRIPTION	DIMENSIONS	
103585	Sterile Liner Set	200 Sterile liners/case	
103604	Dispenser, (wall mounted)	1	
103456	Tray	1	
103457	Tray Lid, clear	1	

"The sterile liner system is one of the most important reasons why many facilities are adopting the CLEANASCOPE™ System – the colour coded liners take the guesswork out of endoscope reprocessing"

Libby Thomson, former Lead Nurse, Bowel Cancer Screening Centre and Nurse Endoscopist

ADVANTAGE PLUSTM and INTERCEPTTM are trademarks of Medivators Inc.

CLEANASCOPE™ is a trademark of Cantel (UK) Limited.

### **MEDIVATORS IS NOW CANTEL**

### www.cantelmedical.eu

Medivators BV | p: +31 45 5 471 471 (customer service)

