
CANTEL REPROCESSING REIMAGINED

IN KEEPING WITH OUR DEDICATION TO INFECTION PREVENTION

Coming Soon: RAPID Workflow — Optimise your reprocessing.

Cross-contamination and instrument damage is reduced with the new integrated workflow of RAPIDAER™ Endoscope Reprocessor, ENDODRY™ Drying and Storage System, and CLEANASCOPE™ ADVANTAGE Transport and Storage System, as endoscope handling is greatly minimised.



RAPIDAER™

Endoscope Reprocessor,
ENDODRY™ Drying
and Storage System
and CLEANASCOPE™
ADVANTAGE Transport
and Short-term Storage
System — the most
advanced and proven
technologies for
increased efficiency and
minimal footprint.



RAPID WORKFLOW

RAPIDAER™ Endoscope Reprocessor

Increase Productivity

- Reprocesses three scopes per hour
- Single bay design for complete install flexibility
- External scope loading for efficient throughput and reduced operator strain
- Compact footprint for maximum room capacity

Guarantee Safety

- Compliant with EN ISO 15883-1, EN ISO 15883-4
- Automated daily Thermal Self Disinfection cycle at >80° C as by EN ISO 15883-4
- Complete traceability through latest RFID technology
- Fully Independent Monitoring System (IMS)
- Remote diagnostics enable fast and efficient technical support and servicing, maximising reprocessor availability



ENDODRY™ Drying and Storage System

Ensure Safety

- Locking doors protect endoscope inventory
- Compliant with EN 16442
- Designed to use air for direct channel drying meeting ISO 8573-1 requirements
- Available in single-sided or pass-through models (front and back door) for separation of wet and dry areas

Increase Efficiency

- Colour touch screen display for quick identification of endoscope storage time
- Barcode scanner for complete traceability
- SMART Connect data interface with RAPIDAER™ Reprocessor and CANEXIS™ Integrated Workflow Solution for increased security
- Dries endoscopes within three hours (internal channels within one hour) and can be stored up to 31 days without an increase in microbial levels



CLEANASCOPE™ ADVANTAGE Transport and Storage System

- Protects your endoscopes by minimising direct handling and reducing contact with airborne contaminants during transport and short-term storage
- Sterile liners provide clear identification of "clean" or "contaminated" scopes
- Robust reusable trays easily fit RAPID Workflow load carrier and help prevent damage to equipment
- Hard lids available for added protection
- Available in 5, 6 and 10 tray options



CONTACT YOUR LOCAL SALES REPRESENTATIVE

All company and product names are trademarks of Cantel Medical Corp., its affiliates or related companies, unless otherwise noted.

www.cantelmedical.eu

CONTACT US FOR MORE INFORMATION

Medivators BV | t: +31 (0)45 547 1471 | e: csdistributors@medivators.com

