



# LiNA Xcise™

*Cordless Morcellator*




**Patent Pending**



# LiNA Xcise™

The LiNA Xcise is a cordless fully disposable morcellator, intended for tissue morcellation during laparoscopic gynaecological procedures. The unique design of the LiNA Xcise provides easy control during hysterectomy and myomectomy procedures.



## safe

### Three-step positioning

- **Secure introduction**

The protective cannula can be set in three positions: "Safety Guard", "Cut 1" and "Cut 2". In position "Safety Guard" the protective cannula serves as a safety guard when introducing the morcellator into the abdomen.

## versatile

### The next generation of morcellation

- **Cordless**

The morcellator can be used right of the packaging, no assembly needed.

- **Self contained**

No capital equipment or additional investment is needed.

- **Mobility**

The fact that no cable is needed for the morcellator facilitates unlimited mobility of the OR staff.

- **Environmental**

LiNA Xcise is RoSH compliant, so the product including battery does not contain any lead or cadmium.



## functional

### Easy morcellation, minimal preparation

- **Three-step positioning**

"Safety Guard" for secure introduction. "Cut 1" for surface peeling. "Cut 2" for deeper coring.

- **Fully hand controlled**

The excellent activation button switches the morcellation on and off instantly.

- **Ergonomic design**

The design of the morcellator complements most hand sizes as well as left- and right-handed surgeons.



### Ordering information:

LiNA Xcise™ **CE 0459**

REF.	Product	Packing
MOR-1515	LiNA Xcise – 15 mm morcellator with obturator – single use	Box of 6, sterile
MOR-1515-1	LiNA Xcise – 15 mm morcellator with obturator – single use	Box of 1, sterile

