

OPTICAL MINI URETHROTOME

E-line

FOR TELESCOPE 3.3 MM, 0°



12 o'clock cutting position

AUTOCLAVE
134° C / 273° F

for optically controlled slitting of urethral strictures.

Product no.

Sheath for urethrotome 16.5 Fr

with supply and drain stopcock,
lateral instrument port (4 Fr),
distal end slanted, and automatic lock



8670.011

Obturator for urethrotome 16.5 Fr

Oval, cannulated, distal end rounded



8670.111

Outer sheath for urethrotome 16.5 Fr

Oval, laterally open, for placing a 16 Fr balloon catheter



8670.061

Outer sheath for urethrotome 16.5 Fr

Oval, laterally open,
for placing a 12 Fr balloon catheter



8670.161

Working element, urethrotome, passive 0°

Receptacle for stricture scalpels



8670.911

Also:



Telescope

Direction of view

Ø

WL

TL

Product no.

Telescope 0° Ø 3.3 mm WL 297 mm

0°

3.3 mm

297 mm

362 mm

8660.424

HD
High Definition

Recommended accessories

Product no.

Perforated basket for standard telescope 2

Internal dimensions (W x H x D): 454 x 52 x 40 mm
External dimensions (W x H x D): 481 x 54 x 59 mm

38016

Fiber light cable bndl

consisting of:
Fiber light cable (80662523), projector-side adapter (8095.07),
and endoscope-side adapter (809509)

Ø 2.5 mm, 2.3 m long

806625231

Stricture scalpels

see URO23

STRICTURE SCALPELS

FOR OPTICAL MINI URETHROTOME 8670

Stricture scalpels with ceramic cutter		Cutting edge	Product no.
 Solid sheath	Scalpel 16.5 Fr ceramic, semi-circular	Semi-circular	8670.961
 Sheath with central bore hole (for 4 Fr ureteral catheter)	Scalpel 16.5 Fr ceramic, semi-circular	Semi-circular	8670.965
 Sheath with central bore hole (for 4 Fr ureteral catheter)	Scalpel 16.5 Fr ceramic, lancet-shaped	Lancet-shaped	8670.981

Recommended accessories		Product no.
Protective sterilization sleeve Receptacle for max. 5 electrodes, reusable	 AUTOCCLAVE 134° C / 273° F	33030