

LiNA Pencils

For open monopolar electrosurgery

LiNA offers a wide range of **Adaptors** for HF generators.



LiNA **Pencil Holster** for safer electrosurgery.



LiNA **Tip Cleaner** for cleaning of electrodes.



Ordering information:

LiNA Pencils

REF:	Product	Packing
	LiNA Surgipen - Reusable - The Surgipen is available with the following cable-lengths: 3 or 5 meters. Indicate when ordering by adding the length to product code, e.g. 508-L-3 for 3 meter cable.	
508-L-3 508-L-5	LiNA Surgipen - Reusable - Autoclavable for 2.35mm electrodes . With 3m silicone cable and standard US-3-pin connector. With 70mm blade electrode.	Box of 10
508-L-3-ACC 508-L-5-ACC	LiNA Surgipen - Reusable - Autoclavable - for 2.35mm electrodes . with silicone cable and ACC/ICC connector for Erbe ACC/ICC generators. With 70mm blade electrode.	Box of 10
512-L-3 512-L-5	LiNA Surgipen - Reusable - Autoclavable - for 4mm hexagonal electrodes . With silicone cable and standard US-3-pin connector.	Box of 10
512-L-3-ACC 512-L-5-ACC	LiNA Surgipen - Reusable - Autoclavable for 4mm hexagonal electrodes . With silicone cable and ACC/ICC connector for Erbe ACC/ICC generators.	Box of 10
619	LiNA Footswitch Pen - foot controlled - Reusable - Autoclavable - for 2.35mm electrodes. With 3m silicone cable and standard 4mm male connector.	Box of 10
516	LiNA Fingerswitch Pen - Single use - for 2.35mm electrodes. With 3m cable and standard US-3-pin connector. With 70mm blade electrode.	Box of 25, sterile
516-H	LiNA Fingerswitch Pen - Single use - for 2.35mm electrodes. With 3m cable and standard US-3-pin connector. With 70mm blade electrode + With Pencil Holster	Box of 25, sterile
517	LiNA Footswitch Pen - foot controlled - Single use - for 2.35mm electrodes. With 3m cable and 4mm standard male connector. With 70mm blade electrode.	Box of 25, sterile
580	LiNA Tip Cleaner - Single use - Self Adhesive. 45x50mm with X-ray detectable thread.	Box of 25, sterile
519-R 519-S	LiNA Pencil Holster - Semi reusable - Autoclavable. LiNA Pencil Holster - Single use.	Box of 10 Box of 25
571-ACC 571-AE 571-B 571-E 571-L 571-M 571-N 571-S1 571-S2 572	LiNA HF generator Adaptor - reusable - for Erbe ACC/ICC generators. LiNA HF generator Adaptor - reusable - for Aesculab generators. LiNA HF generator Adaptor - reusable - for Bovie generators. LiNA HF generator Adaptor - reusable - for Old Erbe generators (Erbotom). LiNA HF generator Adaptor - reusable - for Lamidey generators. LiNA HF generator Adaptor - reusable - for Martin generators. LiNA HF generator Adaptor - reusable - for Neomed generators. LiNA HF generator Adaptor - reusable - for Siemens type 1 generators. LiNA HF generator Adaptor - reusable - for Siemens type 2 generators. LiNA HF generator Adaptor - reusable - 8mm connector for foot controlled pencils. For REF 619 and 517.	Box of 1



REF 508/512



**LiNA Surgipen
Reusable**

With integrated silicone cable. Autoclavable 100 times. Available with standard US-3-pin connector or ACC/ICC connector. Available for 2.35mm electrodes and for 4mm Hexagonal electrodes. With 3 or 5m cable.

REF 619



**LiNA Footswitch Pen
Reusable**


Foot controlled pencil with integrated silicone cable. Autoclavable. Available for 2.35mm electrodes and with 3m cable with 4mm male connector. Adaptor with 8mm male connector available REF-572.

LiNA Pencils

For monopolar electrosurgery procedures

LiNA offers a range of reusable and single use Pencils for monopolar electrosurgery. Both finger- and footswitch controlled electrosurgical pencils are available with integrated cables. All LiNA Pencils are cost effective and offers optimal comfort and manoeuvrability for the surgeon.

REF 516



**LiNA Fingerswitch Pen
Single Use**

With 3m cable and US-3-pin connector. Available for 2.35mm electrodes.

REF 517



**LiNA Footswitch Pen
Single Use**

Foot controlled pencil with integrated cable. Available for 2.35mm electrodes and with 3m cable with 4mm male connector. Adaptor with 8mm male connector available REF-572.