



OverStitchTM Endoscopic Suturing System

Surgical Suturing on
Your Flexible Endoscope,
Less Invasive Options
for Your Patient.



Defining the Future
of Endoscopic SurgeryTM

OverStitch™ Endoscopic Suturing System



Full Thickness Surgical Suturing Delivered Through a Flexible Endoscope

- Capable of deploying multiple running and interrupted stitches with a single insertion of the endoscope
- Unique design mimics curved needle, which allows for controlled depth of suture placement
- Available in both absorbable and non-absorbable suture materials:
 - 2-0 Polypropylene (PP)
 - 2-0 Polydioxanone (PDO)
- Mounts quickly and easily on a therapeutic endoscope and allows sutures to be deployed under direct visualization
- Single-handed operation allows physician to maintain system control
- Reloadable while maintaining visibility of the operative site
- Unique cinch design provides secure closure without the need to tie surgical knots
- Single-patient use for maximum reliability and patient safety



Product Family	Code	Description	Package Qty
OverStitch™	ESS-G02-160	OverStitch™ Endoscopic Suturing System	3 units
	PLY-G02-020	OverStitch™ 2-0 Polypropylene Suture	12 units
	PDO-G02-020	OverStitch™ 2-0 Polydioxanone Suture	12 units
	CNH-G01-000	OverStitch™ Suture Cinch	6 units
Accessories	THX-165-028	Tissue Helix, 165cm Length, 2.8mm OD	6 units
	OVT-027-160	OverTube™ Endoscopic Access System	3 units

Defining the Future
of Endoscopic Surgery™

www.apolloendo.com



Apollo Endosurgery, Inc. • 7000 Bee Caves Road, Suite 350 • Austin, TX 78746 • Tel. 877-ENDO-130 / 512-328-9990 • www.apolloendo.com

OverStitch™ and Defining the Future of Endoscopic Surgery™ are trademarks of Apollo Endosurgery, Inc. Copyright ©2011 Apollo Endosurgery, Inc. All rights reserved.
Caution: Federal law (USA) restricts this device to sale by or on the order of a physician. MKT-00015-00R04/060112