





Daisy

The nitinol Daisy basket is designed for use in the kidney and ureter, providing access and predictable stone capture.

The Daisy basket is available in 4 and 6 wire.



Stonefree
Solutions
Evolution in urology

About us

With a focus on delivering quality proven clinical solutions and developing value added innovation, our entire team has vast experience in the Urology stone management field. We have the 'know how', that gives you the confidence to choose Stone Free for your patients.



Daisy Ordering Information

Product Name	Part Number	Fr	Basket Diameter (mm)	Basket Configuration	Length (cm)
Tipless Stone Extractor	SFS1500DAI1	1.5	10	4 wire, tipless	120
	SFS1500DAI2	1.9	10	4 wire, tipless	120
	SFS1500DAI3	2.2	10	4 wire, tipless	120
	SFS1500DAI4	2.2	20	4 wire, tipless	120
	SFS1500DAI5	3	10	4 wire, tipless	120
	SFS1500DAI6	3	20	4 wire, tipless	120
	SFS1500DAI8	4.5	20	4 wire, tipless	120
	SFS1500DAI7	1.9	10	6 wire, tipless	120