





Tulip



Designed to provide easy and controlled positioning of stones in the kidney and ureter.

The Tulip basket's nitinol hand formed based design provides safe, kink resistant and predictable performance.



Stonefree
Solutions
Evolution in urology

About us

With a focus on delivering quality proven clinical solutions and developing value added innovation, our entire team has vast experience in the Urology stone management field. We have the 'know how', that gives you the confidence to choose Stone Free for your patients.



Tulip Ordering Information

Product Name	Part Number	Fr	Basket Diameter (mm)	Basket Configuration	Length (cm)
3 Wire Stone Extractor	SFS1900TUL1	1.7	8	3 wire	120
	SFS1900TUL2	1.7	11	3 wire	120
	SFS1900TUL3	1.9	10	3 wire	120
	SFS1900TUL4	2.2	8	3 wire	120
	SFS1900TUL5	2.2	11	3 wire	120