

HAWK

Mobile Diagnostics
and Treatment



otoprnt®

HAWK

The smallest mobile ENT treatment unit.

Developed for flexible use in bedside visits, outpatient clinics, nursing homes and armee field container. With it's stable and light weight aluminium casing the unit can be easily transported in the car trunk. The instruments can be stored save in the lockable storage doors. The HAWK combines the functional comfort of a full equipped ENT unit with smalest space requirements and mobility.

Technical Data

Suction unit:	24l / min, 90% max vacuum, integrated air filter, vacuum gauge, double overflow fuse
Halogen light source:	1 x 100 W halogen outlet with two layer thermal shield
Quivers:	2 quivers for rigid endoscopes, 1 quiver for flexible endoscopes, head mirror rest
Compressed air unit:	0.5-3.0 bar adjustable, with fine regulation handpiece, 2 spray bottles
Ear rinsing device:	38°C water temperature, autarchic system with pressure pump, including autoclavable filter unit with filter cartridge, detachable 1.5l water tank with cap
Hygiene:	all functional parts with contact to the patient are autoclavable
Casing:	light aluminium casing, lacquered in white color, other colours on request
Mobility:	4 antistatic wheels, 2 guide wheels in front, 4 handles for transport, auto furl power cable
Weight:	32.5 kg
Power supply:	110 - 230V, 50 - 60 Hz

All designs and specifications are subject to change without notice.



Happersberger otopront GmbH
Langgasse 90
D-65329 Hohenstein

Tel. +49 (0) 61 20 - 92 17 0
Fax. +49 (0) 61 20 - 92 17 60

www.otopront.de
E-Mail: info@otopront.de